

# Technical Information

## Rxn-10 Raman spectroscopic probe

A versatile probe for your Raman spectroscopy needs

### Application

Designed for product and process development, the Rxn-10 probe is trusted to deliver high-performance measurements across a wide spectral range. Compact, lightweight, and flexible, it is ideal for both solids and liquids analysis. With interchangeable optics, it adapts easily to diverse applications. Now compatible with our new KFOC1B Raman fiber-optic cable, it offers enhanced certifications and greater installation flexibility for lab and industrial environments.

- **Chemical:** reaction monitoring, blending, catalyst monitoring, hydrocarbon speciation, process unit optimization.
- **Polymer:** polymerization reaction monitoring, extrusion monitoring, polymer blending
- **Pharmaceutical:** active pharmaceutical ingredient (API) reaction monitoring, crystallization
- **Biopharmaceutical:** cell culture and fermentation monitoring, optimization, control
- **Food and beverage:** zonal heterogeneity mapping of meats and fish

### Device properties

- 6061 aluminum, 316L stainless steel, and 303 stainless steel
- PVC jacketed, proprietary construction
- proprietary electro-optical (EO), or FC to EO fiber converter(s) for non-embedded systems

### Your benefits

- Multipurpose use for both solids and liquids measurement
- Lightweight and compact
- Integrated laser safety interlock, including "laser on" indication and probe shutter
- Flexible output compatible with a range of sampling options
- Easy switching of non-contact, immersion, and bioprocessing optics to suit a variety of applications
- Wide spectral range, including access to the critical low-wavenumber region
- Upgraded, CMR-certified KFOC1B Raman fiber-optic cable option for improved fire resistance, simplified regulatory compliance, and enhanced flexibility for easier routing and handling







# Table of Contents

<b>About this document .....</b>	<b>4</b>	Installation .....	7
Symbols .....	4	<b>Specifications .....</b>	<b>8</b>
<b>Function and system design .....</b>	<b>5</b>	Probe specifications.....	8
Application.....	5	Fiber-optic cable specifications .....	9
Laser safety interlock.....	5	Probe dimensions .....	10
Rxn-10 probe .....	5	MPE: ocular exposure.....	11
Rxn-10 probe optics .....	6	MPE: skin exposure .....	11

## About this document

### Symbols

### Safety symbols

<p> <b>WARNING</b></p> <p><b>Causes (/consequences)</b> Consequences of noncompliance (if applicable) ► Corrective action</p>	<p><b>Standard precautions for laser products should be observed.</b></p> <ul style="list-style-type: none"> <li>► Probes should always be shuttered or pointed away from people toward a diffuse target if not installed in a sample chamber.</li> </ul>
<p> <b>CAUTION</b></p> <p><b>Causes (/consequences)</b> Consequences of noncompliance (if applicable) ► Corrective action</p>	<p><b>The laser input into the Rxn-10 probe must not exceed 499 mW.</b></p> <p><b>If stray light is allowed to enter an unused probe, it will interfere with data collected from a used probe and may cause calibration failure or measurement errors.</b></p> <ul style="list-style-type: none"> <li>► Unused probes should ALWAYS be shuttered to prevent stray light from entering the probe. If an optic cap is available, place it on the unused optic.</li> </ul>
<p><b>NOTICE</b></p> <p><b>Cause/situation</b> Consequences of noncompliance (if applicable) ► Action/note</p>	<p><b>When installing the probe head <i>in situ</i>, the user must ensure that there is strain relief at the installation location which complies with bend radius specifications.</b></p>

## Function and system design

**Application** Use of the device for any purpose other than that described, poses a threat to the safety of people and of the entire measuring system and invalidates any warranty.

**Laser safety interlock** The Rxn-10 probe, as installed, forms part of the interlock circuit. If the fiber-optic cable is severed, the laser will turn off within milliseconds of the breakage.

**NOTICE**

**Permanent damage may result if cables are not routed appropriately.**

- ▶ Handle probes and cables with care, ensuring they are not kinked.
- ▶ Install fiber-optic cables with a minimum bend radius according to the *Raman fiber-optic cable Technical Information (TI01641C)*.

**Rxn-10 probe**










Figure 1: Rxn-10 probe

#	Name	Description
1	Fiber-optic cable	Connects the probe to the Raman Rxn analyzer via the electro-optical (EO) fiber-optic cable attached to the Rxn-10 probe.
2	Laser emission indicator	When there is potential for the laser to be energized, the indicator light is illuminated.
3	Laser beam shutter	Can be closed to prevent laser emission. Position "I" indicates emission potential. Moving the lever past position "O" indicates emission is shuttered.
4	Thumb screw	Tighten to secure optics onto the probe when there is not a threaded interface.
5	Optics interface	Insert optics or threaded adapter.

**Rxn-10 probe optics**

The probe is compatible with the following optics to meet the requirements of different applications:

Optics		Applications
Non-contact optics	 <p>A0048410 A0048676</p>	For use with solids or turbid media. Also well-suited for delicate or corrosive liquids when sample contamination or damage to optical components is a concern.
Immersion optics (IO)	 <p>A0048411</p>	For use in reaction vessels, laboratory reactors, or process streams.
bio-Optic	 <p>A0048412</p>	For use with continuous inline measurement in benchtop bioreactor/fermenter applications requiring headplate entry.
Bio multi optic and bio sleeve	 <p>A0051184</p>	For use with continuous inline measurement in benchtop bioreactor/fermenter applications requiring headplate entry.
Raman optic system for single use	 <p>A0048413</p>	For use with disposable fittings for single-use applications.

Optics		Applications
Raman flow assembly (includes micro flow bench and micro flow cell)		For use with lower flow rate liquids, where monitoring a dynamic process stream provides valuable information, and speed or limit of detection are particularly important.

Installation

During installation, standard eye and skin safety precautions for Class 3B laser products (as per EN-60825/IEC 60825-14 or ANSI Z136.1) should be observed.

## Specifications

### Probe specifications

Specifications for the Rxn-10 probe are listed below.

Item		Description
Laser wavelength	with non-contact or immersion optic	532 nm, 785 nm, or 1000 nm
	with bIO-Optic or Raman optic system for single use	785 nm or 1000 nm
	with bio multi optic and bio sleeve or micro flow bench and micro flow cell	785 nm
Maximum laser power into probe head		< 499 mW
Working distance		Refer to <i>Accessory optics for the Rxn-10 probe Technical Information (TI01635C)</i>
Sample interface		Refer to <i>Accessory optics for the Rxn-10 probe Technical Information (TI01635C)</i>
Polarization at sample		Unpolarized
Ambient temperature		−10 to 70 °C (14 to 158 °F)
Temperature ramp		≤ 30 °C/min (≤ 54 °F/min)
Relative humidity		20 to 60 %, non-condensing
Spectral coverage		Probe spectral coverage is limited by the coverage of the analyzer being used
Laser power at sample	532 nm (with standard 120-mW laser)	> 45 mW
	785 nm (with standard 400-mW laser)	> 150 mW
	1000 nm (with standard 400-mW laser)	> 150 mW
Materials of construction	probe body	6061 aluminum, 316L stainless steel, and 303 stainless steel
	fiber-optic cable	Design: PVC jacketed, proprietary construction Connections: proprietary electro-optical or FC to EO fiber converter(s) for non-embedded systems
Probe	length (not including fiber-optic cable bend radius)	203 mm (8 in)
	length (including fiber-optic cable bend radius)	356 mm (14.02 in)
	diameter (not including cable)	19 mm (0.75 in)
	weight (including cable)	0.5 kg (approximately 1 lb)



**Fiber-optic cable specifications** Specifications for the fiber-optic cables are listed below.

<b>KFOC1B Raman fiber-optic cable</b>	
<b>Item</b>	<b>Description</b>
General features	Integrated copper conductor wire for interlock capability Aramid (Kevlar) internal strength members Flame retardant Fungus resistant
Cable rating (cable only)	Operating temperature: $-40^{\circ}\text{C}$ to $70^{\circ}\text{C}$ ( $-40^{\circ}\text{F}$ to $158^{\circ}\text{F}$ ) Storage temperature: $-55^{\circ}\text{C}$ to $70^{\circ}\text{C}$ ( $-67^{\circ}\text{F}$ to $158^{\circ}\text{F}$ ) Certified: CSA-C/US AWM I/II, A/B, 80C, 30V, FTI, FT2, VW-1, FT4 Rated: AWM I/II A/B 80C 30V FT4
Bend radius	152.4 mm (6 in)
Termination	Electro-optical (EO) with connectors


The KFOC1B Raman fiber-optic cable features an improved rating, CMR-certification, ensuring easier compliance with local laws and regulations. This certification supports smoother implementation in process environments. Independently tested and certified by a third party, these cables offer enhanced protection against the spread of fire.

With the CMR rating, the KFOC1B Raman fiber-optic cable is ready for immediate installation in cable trays, risers, and all conduit types with no additional assessments required.

<b>KFOC1B Raman fiber-optic cable</b>	
<b>Item</b>	<b>Description</b>
General features	Integrated copper conductor wire for interlock capability Fiber-reinforced plastic (FRP) strength members Flame retardant Fungus resistant
Cable rating (cable only)	Operating temperature: $-40^{\circ}\text{C}$ to $70^{\circ}\text{C}$ ( $-40^{\circ}\text{F}$ to $158^{\circ}\text{F}$ ) Storage temperature: $-55^{\circ}\text{C}$ to $70^{\circ}\text{C}$ ( $-67^{\circ}\text{F}$ to $158^{\circ}\text{F}$ ) Certified: cULus AWM I/II, A/B, 80C, 30V, FTI, FT2, VW-1, FT4 Rated: CMR-FO, AWM I/II A/B 80C 30V FT4
Bend radius	152.4 mm (6 in)
Termination	Electro-optical (EO) connectors

**Probe dimensions**

The dimensions for the Rxn-10 probe are shown below.



A0048400

Figure 2. Rxn-10 probe dimensions

Dimension	Measurement	Description
L1	111 mm 4.37 in	Length of probe body without fiber-optic cable
L2	203 mm 8 in	Length with fiber-optic cable connected Note: This does not include additional minimum bend radius of cable
H	33 mm 1.3 in	Height of probe including thumb screw
Ø	19 mm 0.75 in	Diameter of probe, not including fiber-optic cable

### MPE: ocular exposure

Refer to the tables below from the ANSI Z136.1 standard to calculate the maximum permissible exposure (MPE) for point source ocular exposure to a laser beam.

A correction factor ( $C_A$ ) may also be required and can be determined below.

Wavelength $\lambda$ (nm)	Correction factor $C_A$
400 to 700	1
700 to 1050	$10^{0.002(\lambda-700)}$
1050 to 1400	5

Maximum permissible exposure (MPE) for point source ocular exposure to a laser beam			
Wavelength $\lambda$ (nm)	Exposure duration $t$ (s)	MPE calculation	
		(J·cm <sup>-2</sup> )	(W·cm <sup>-2</sup> )
532	$10^{-13}$ to $10^{-11}$	$1.0 \times 10^{-7}$	-
	$10^{-11}$ to $5 \times 10^{-6}$	$2.0 \times 10^{-7}$	-
	$5 \times 10^{-6}$ to 10	$1.8 t^{0.75} \times 10^{-3}$	-
	10 to 30,000	-	$1 \times 10^{-3}$

Maximum permissible exposure (MPE) for point source ocular exposure to a laser beam				
Wavelength $\lambda$ (nm)	Exposure duration $t$ (s)	MPE calculation		MPE where $C_A = 1.4791$
		(J·cm <sup>-2</sup> )	(W·cm <sup>-2</sup> )	
785 and 993	$10^{-13}$ to $10^{-11}$	$1.5 C_A \times 10^{-8}$	-	$2.2 \times 10^{-8}$ (J·cm <sup>-2</sup> )
	$10^{-11}$ to $10^{-9}$	$2.7 C_A t^{0.75}$	-	Insert time (t) and calculate
	$10^{-9}$ to $18 \times 10^{-6}$	$5.0 C_A \times 10^{-7}$	-	$7.40 \times 10^{-7}$ (J·cm <sup>-2</sup> )
	$18 \times 10^{-6}$ to 10	$1.8 C_A t^{0.75} \times 10^{-3}$	-	Insert time (t) and calculate
	10 to $3 \times 10^4$	-	$C_A \times 10^{-3}$	$1.4971 \times 10^{-3}$ (W·cm <sup>-2</sup> )

### MPE: skin exposure

Refer to the table below from the ANSI Z136.1 standard to calculate the MPE for skin exposure to a laser beam.

Maximum permissible exposure (MPE) for skin exposure to a laser beam				
Wavelength $\lambda$ (nm)	Exposure duration $t$ (s)	MPE calculation		MPE where $C_A = 1.4791$
		(J·cm <sup>-2</sup> )	(W·cm <sup>-2</sup> )	
532, 785 and 993	$10^{-9}$ to $10^{-7}$	$2 C_A \times 10^{-2}$	-	$2.9582 \times 10^{-2}$ (J·cm <sup>-2</sup> )
	$10^{-7}$ to 10	$1.1 C_A t^{0.25}$	-	Insert time (t) and calculate
	10 to $3 \times 10^4$	-	$0.2 C_A$	$2.9582 \times 10^{-1}$ (W·cm <sup>-2</sup> )

[www.addresses.endress.com](http://www.addresses.endress.com)

---