

# Technical Information

## Rxn-46 Raman spectroscopic probe

### System design and specifications

#### Application

The union between our Raman analyzers equipped with Rxn-46 bioprocess probe technology and the BioPAT® Spectro platform by Sartorius offers the market an ideal interface to high-throughput development through single-use commercial manufacturing.

Recommended cell culture applications include glucose, lactate, amino acids, cell density, titer, and more.

#### Device properties

Our Raman bioprocessing probe technology has been adapted to fit Sartorius's BioPAT® Spectro platform, utilizing the same probe design for Ambr® 15, Ambr® 250, and Biostat STR® bioreactors.

#### Your benefits

- Enables faster, easier, and more robust model building via integration with Ambr® 15 and Ambr® 250
- Allows high-throughput process development which supports Quality by Design (QbD)
- Provides a more efficient transfer to Biostat STR® for single-use manufacturing
- Offers a scale-independent interface from 15 mL in the laboratory to 2000 L in the production suite
- Requires no probe cleaning, sterilization, or frequent maintenance due to non-contact sampling



## Table of Contents

<b>Function and system design .....</b>	<b>3</b>
Application .....	3
Rxn-46 probe.....	3
Laser safety interlock .....	4
Installation .....	4
Analyzer compatibility .....	5
<b>Specifications .....</b>	<b>6</b>
General specifications .....	6
Probe dimensions: side view.....	6
Probe dimensions: top view .....	7
MPE: ocular exposure.....	7
MPE: skin exposure .....	7

## Function and system design

### Application

Using the device for any purpose other than that described may compromise personal safety, damage the measuring system, and invalidates any warranty.

### Rxn-46 probe




Figure 1. Rxn-46 probe

#	Description
1	Connection to process equipment
2	Probe slider in Operate position
3	Probe body
4	Spring-loaded fiber connector cap
5	Fiber cable
6	Fiber cable connector
7	Fiber cable connector clip

## Laser safety interlock

The Rxn-46 probe, as installed, forms part of the interlock circuit. The interlock circuit is a low-current electrical loop. If the fiber cable is severed, the laser will turn off within milliseconds of the breakage.

### NOTICE

**Permanent damage may result if cables are not routed appropriately.**

- ▶ Handle probes and cables with care, ensuring they are not kinked.
- ▶ Install fiber cables with a minimum bend radius according to the *Raman fiber-optic cable Technical Information (TI01641C)*.

The interlock connector in the fiber cable must be plugged into the interlock socket on a Raman Rxn analyzer, and is automatically connected when the fiber-optic cable process connector is plugged into the Rxn-46 probe. When there is potential for the laser to be energized, the laser interlock indicator light on the probe body is illuminated.




Figure 2. Location of laser interlock indicator light (1)

## Installation

The Rxn-46 probe only interfaces to Sartorius's BioPAT® Spectro compatible parts.

During installation, standard eye and skin safety precautions for Class 3B laser products (as per EN 60825/IEC 60825-14) should be observed. Additionally, observe the following:

<b>⚠ WARNING</b>	<p><b>Standard precautions for laser products should be observed.</b></p> <p>When not installed in a sample chamber, probes should always be capped, pointed away from people, and pointed toward a diffuse target.</p>
<b>⚠ CAUTION</b>	<p><b>If stray light is allowed to enter an unused probe, it will interfere with data collected from a probe in-use and may cause calibration failure or measurement errors.</b></p> <p>Unused probes should ALWAYS be capped to prevent stray light from entering the probe.</p>
<b>NOTICE</b>	<p><b>When installing the probe <i>in situ</i>, the user must provide strain relief to the fiber-optic cable at the probe installation location.</b></p>

**Analyzer compatibility**

The Rxn-46 probe is compatible with the Endress+Hauser Raman Rxn analyzers below operating at 785 nm.

- Ambr<sup>®</sup> 15 and Ambr<sup>®</sup> 250: Raman Rxn2 analyzer, single-channel, benchtop
- Biostat STR<sup>®</sup>: Raman Rxn2 or Rxn4 analyzers, up to four channels; benchtop or mobile wheeled cart (Raman Rxn2); rack mounted or NEMA 4x enclosure (Raman Rxn4)

## Specifications

### General specifications

General specifications for the Rxn-46 probe are listed below.

Item	Description
Laser wavelength	785 nm
Spectral coverage	probe spectral coverage is limited by the coverage of the analyzer being used
Maximum laser power into probe	< 499 mW
IEC 60529 for (EO) right angle connector	IP65
North American TYPE rating for (EO) right angle connector	TYPE 13 <sup>1</sup>
Probe operating temperature	10 to 50 °C (probe is non-contact) (50 to 122 °F)
Probe dimensions (standard)	162 x 159 x 52 mm (6.4 x 6.3 x 2.0 in)

<sup>1</sup> This is a self-declaration of conformance to UL 50E TYPE 13 requirements. It does not constitute UL certification or authorization to use the UL mark. The Rxn-46 Raman spectroscopy probe optics are not air or watertight, therefore we do not make any environmental rating claims on that portion of the probe.

All fiber-optic cable specifications can be found in the *Raman fiber-optic cables KFOC1 and KFOC1B Technical Information (TI01641C)*.

### Probe dimensions: side view




Figure 3. Rxn-46 probe side view. Dimensions: mm (in)

**Probe dimensions: top view**




Figure 4. Rxn-46 probe top view. Dimensions: mm (in)

**MPE: ocular exposure**

Refer to the tables below from the ANSI Z136.1 standard to calculate the maximum permissible exposure (MPE) for point source ocular exposure to a laser beam.

A correction factor ( $C_A$ ) may also be required and can be determined below.

Wavelength $\lambda$ (nm)	Correction factor $C_A$
400 to 700	1
700 to 1050	$10^{0.002(\lambda-700)}$
1050 to 1400	5

MPE for point source ocular exposure to a laser beam				
Wavelength $\lambda$ (nm)	Exposure duration $t$ (s)	MPE calculation		MPE where $C_A = 1.4791$
		(J·cm <sup>-2</sup> )	(W·cm <sup>-2</sup> )	
785	$10^{-13}$ to $10^{-11}$	$1.5 C_A \times 10^{-8}$	-	$2.2 \times 10^{-8}$ (J·cm <sup>-2</sup> )
	$10^{-11}$ to $10^{-9}$	$2.7 C_A t^{0.75}$	-	Insert time (t) and calculate
	$10^{-9}$ to $18 \times 10^{-6}$	$5.0 C_A \times 10^{-7}$	-	$7.40 \times 10^{-7}$ (J·cm <sup>-2</sup> )
	$18 \times 10^{-6}$ to 10	$1.8 C_A t^{0.75} \times 10^{-3}$	-	Insert time (t) and calculate
	10 to $3 \times 10^4$	-	$C_A \times 10^{-3}$	$1.4971 \times 10^{-3}$ (W·cm <sup>-2</sup> )

**MPE: skin exposure**

Refer to the table below from the ANSI Z136.1 standard to calculate the MPE for skin exposure to a laser beam.

MPE for skin exposure to a laser beam				
Wavelength $\lambda$ (nm)	Exposure duration $t$ (s)	MPE calculation		MPE where $C_A = 1.4791$
		(J·cm <sup>-2</sup> )	(W·cm <sup>-2</sup> )	
785	$10^{-9}$ to $10^{-7}$	$2 C_A \times 10^{-2}$	-	$2.9582 \times 10^{-2}$ (J·cm <sup>-2</sup> )
	$10^{-7}$ to 10	$1.1 C_A t^{0.25}$	-	Insert time (t) and calculate
	10 to $3 \times 10^4$	-	$0.2 C_A$	$2.9582 \times 10^{-1}$ (W·cm <sup>-2</sup> )

[www.addresses.endress.com](http://www.addresses.endress.com)

---