Anleitung / Instruction

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## **Kits Filter and Measuring Detector OUSAF44/46**

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**This instruction applies to following spare part kits:**

71108537	Kit OUSAF44 Interference Filter 254 nm
71108539	Kit OUSAF44 Interference Filter 280 nm Re
71108538	Kit OUSAF44 Interference Filter 280 nm Mea
71108536	Kit OUSAF44 Interference Filter 295 nm
71108540	Kit OUSAF44 Interference Filter 302 nm
71108542	Kit OUSAF44 Interference Filter 365 nm
71108541	Kit OUSAF44 Interference Filter 313 nm
71108550	Kit OUSAF44 Measuring Detector Assembly

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## Contents

	Page
1. General notes and notes on safety	2
2. Important note on EasyCal™ systems	3
3. Note on filter usage	4
4. Sensor design / sectional drawing	4
5. Scope of supply	5
6. Replacement instructions	
6.1 Required tools and material	6
6.2 Measurement filter replacement	6
6.3 Measuring detector replacement	8
6.4 Reference filter replacement	9

## 1. General notes

**Please pay attention to the following warnings for your own safety! Please pay attention to the operating instructions of your OUSAF44!**

### 1.1 Qualification

The replacement of the OUSAF4x interference filters and measuring detector may only be done by qualified and authorized personnel in compliance with these instructions.

### 1.2 Usage of the kit

The parts of this service kit are only to be used as spare parts for sensors of the OUSAF44 / OUSAF46 sensors family. Any other use is not permitted.

### 1.3 Electric safety

The lamp of the OUSAF4x sensor operates internally with high voltage. In order to avoid electric shock all cables and power should be disconnected from the sensor prior to performing any maintenance or repair that requires the lamp housing to be opened.

### 1.4 Eyes protection

The OUSAF44 lamp produces UV light. Whenever maintenance is performed on the sensor, remove all cables and power to the sensor. Never view the lamp directly without proper eye protection when it is powered on!

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### 1.5 Burn risk

The lamp heats up during normal operation, so allow 30 minutes for the lamp to cool down after powering it down. Do not touch the lamp with your fingers.

### 1.6 Safety rules

The sensor must be disconnected from the transmitter before starting the work. Read the complete instructions, and be sure to completely understand them before starting the work.

### 1.7 Feedback on the process

Consider the feedback to your process before switching off an active instrument! Observe your process if it is controlled with the help of the involved measuring instrument via switching contacts, via analog signal outputs, or via a field bus communication interface.

### 1.8 Removal of sensors out of the process

For your own safety, pay attention to process pressure, process temperature as well as aggressiveness or toxicity of the medium, if a sensor has to be removed from the process for maintenance or service.

## 2. Important note regarding EasyCal™ systems

If your OUSAF44 Sensor is equipped with an EasyCal™ calibration tool, its certification will be voided if the measuring filter, detector, or any other element of the detector assembly are replaced or otherwise handled.

If your EasyCal™ detector assembly needs to be serviced, you can either replace it using a complete EasyCal™ retrofit kit or contact your local Endress+Hauser service organization to arrange for its repair or recertification.

This instruction is meaningless in this case.

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### 3. Notes on usage of the interference filters


- The filters from the "Interference Filter" kits can be used both as measurement filter or as reference filters. Measuring filter replacement see chapter 6.2, reference filter replacement see chapter 6.4.
- Important:**  
The filter for 280 nm is an exception. Here, the interference filter 280 nm Re (Kit 1108539) is used as reference filter. As measuring filter, the interference filter 280 nm Mea is used (Kit 1108538).
- Observe the installation direction of the filters in chapters 6..2 and 6.4!

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### 4. Sensor design/sectional drawing

The sectional drawing (figure 1) shows the internal design of the OUSAF44 sensor.

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## 5. Scope of supply

### 71108550 Kit OUSAF44 Measuring Detector Assembly

This kits contains the following parts:

- 1 pc. Measurement detector with connector included for OUSAF44 sensors
- 3 pcs. Slotted screws M3-08
- 4 pcs. Split washers M3
- 1 pc. Connector cover
- 1 pc. Instruction for service kit

### 71108538 KIT OUSAF44 Interference Filter 280nm Mea

This kits contains the following parts:

- 1 pc. 280 nm Interference filter Mea (measurement filter)
- 1 pc. O-ring 2-019 Silicone for the detector assembly housing
- 1 pc. Instruction for service kit

### 71108539 KIT OUSAF44 Interference Filter 280nm Re

This kits contains the following parts:

- 1 pc. 280 nm Interference filter Re (reference filter)
- 1 pc. O-ring 2-025 Black Buna (lamp housing)
- 1 pc. O-ring 2-026 Black Buna (lamp adapter)
- 1 pc. O-ring 2-012 Black Buna (reference filter)
- 1 pc. Instruction for service kit

### All other filter kits listed

- 1 pc. Interference filter for the wavelength stated in the kit description
- 1 pc. O-ring 2-025 Black Buna (lamp housing)
- 1 pc. O-ring 2-026 Black Buna (lamp adapter)
- 1 pc. O-ring 2-012 Black Buna (reference filter)
- 1 pc. O-ring 2-019 Silicon (not used for ATEX/FM lamps)
- 1 pc. Instruction for service kit

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## 6. Replacement instructions

### 6.1 Required tools and material


1. 1.5mm Allen key
2. 3mm slotted screwdriver
3. No. 1 Philips Screwdriver
4. 2 adjustable pipe or strap wrenches with minimum opening from 0" to 3"
5. 1 Adjustable wrench with minimum opening from 0" to 3"
6. Talc-free Latex gloves resp. cotton gloves
7. Kit OUA260 lamp and detector O-ring, order number 71142537

### 6.2 Measurement filter replacement

1. Remove the lamp and detector assemblies from the flowcell (1 and 2 in figure 2) by turning them counterclockwise. If needed, use the strap or pipe wrenches, making sure to protect the surface of the sensor in order to prevent deep scratches on it.
2. Remove the top cover of the detector assembly. You may have to use the pipe or strap wrench, as well as the adjustable wrench, to remove it by holding the top cover by the connector (10 in figure 3) with the adjustable wrench and using one of the strap or pipe wrenches to turn the lamp's assembly body. Using the wrenches to open the lamp assembly, please make sure to protect the surface of the sensor in order to prevent deep scratches on it.
3. Using the 1.5mm Allen key, loosen the measurement filter set-screw (8 in figure 3).
4. Replace the measurement filter (9 in fig. 3). Please note the orientation of the filter so you can install the new filter the same manner. Do not touch the filter with bare hands!

#### NOTICE

The arrow on the filter represents the light path. The side with the smaller aperture faces the light source.

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- Reassemble the sensor by following steps 1 to 4 in reverse order and replacing the detector housing O-ring (2-019 Silicone, 5 in figure 3), as well as the lamp and detector assembly O-rings (2-022, 3 in figure 2) with new O-rings. Make sure the lamp assembly is mounted on the side of the flowcell with the shorter window.
- Follow your transmitter procedure to calibrate the measuring loop.

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### 6.3 Measurement detector replacement

1. Remove the lamp and detector assemblies from the flowcell (1 and 2 in figure 2) by turning them counterclockwise. If needed, use the strap or pipe wrenches, making sure to protect the surface of the sensor in order to prevent deep scratches on it.
2. Remove the 4 screws that keep the connector of the measuring detector in place (6 in figure 3).
3. Pull out the connector (7 in figure 3). The measuring detector is attached to it and will come out with it.
4. Reassemble the sensor by following steps 1 to 4 in reverse order and replacing the measuring detector and the lamp and detector assembly O-rings (2-022, 3 in figure 2) with new parts.
5. Follow your transmitter procedure to calibrate the measuring loop.

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Figure 3:  
Detector assembly components

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#### 6.4 Reference filter replacement

1. Remove the lamp and detector assemblies from the flowcell (1 and 2 in figure 2) by turning them counterclockwise. If needed, use the strap or pipe wrenches, making sure to protect the surface of the sensor in order to prevent deep scratches on it.
2. Remove the top cover of the lamp assembly. You may have to use the pipe or strap wrench, as well as the adjustable wrench to remove it by holding the top cover by the connector (11 in figure 4) with the adjustable wrench and using one of the strap or pipe wrenches to turn the lamp's assembly body. If you have to use the wrenches to open the lamp assembly, please make sure to protect the surface of the sensor in order to prevent deep scratches on it.
3. Locate the lamp's holder and reflector (12 in figure 4) and, on it, the reference detector board (17 in figure 5).

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4. Using a slotted screwdriver, remove the screws that hold the reference detector board in place (16 in figure 5).
5. Remove the O-ring of the reference filter (14 in figure 5) and the reference filter (15 in figure 5) by turning the lamp assembly so they can fall out. Please note the orientation of the filter so you can install the new filter the same way.

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**NOTICE**

Do not touch the reference filter with your bare hands!

The arrow on the filter represents the light path. The side with the smaller aperture faces the light source. Insert the filter "arrow showing upwards" (15 in figure 5).

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Figure 5:  
Reference detector and filter

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[https://portal.endress.com/picpool/sft/advicepdfs/71210425ousaf44\\_filter\\_detector.pdf](https://portal.endress.com/picpool/sft/advicepdfs/71210425ousaf44_filter_detector.pdf)

If you need more information about your sensor OUSAF44:

[www.endress.com/OUSAF44](http://www.endress.com/OUSAF44)

Here you'll find (as far as applicable or available for this product):

- Product specifications
- Technical Information (TI)
- Operating instruction (BA)
- Short instruction (KA)
- Certificates
- Accessories
- Maintenance parts
- Spare parts
- and more

Here you will get the original edition of the present instruction: