Technical Information **Memosens CPL57E**

pH measurement for laboratory measurements and random sampling



Digital with Memosens 2.0 technology pH sensor for pure and ultrapure water

Application

Measurements in low-conductivity media, such as boiler feedwater, pure and ultrapure water

Your benefits

- Memosens gel compact pH sensor with easy-to-clean glass body
- Reference system with salt packing for drift-free measurement in low-conductivity media
- Integrated NTC 30K temperature sensor for effective temperature compensation
- Suitable for use with Liquiline Mobile, Liquiline To Go and Memobase Plus

Other advantages of Memosens technology

- Maximum analysis safety with non-contact, inductive signal transmission
- Data security thanks to digital data transmission
- Very easy to use as sensor data saved in the sensor



Function and system design


Measuring principle

pH measurement

The pH value is used as a unit of measurement for the acidity or alkalinity of a medium. The membrane glass of the electrode delivers an electrochemical potential that depends on the pH value of the medium. This potential is generated by the selective accumulation of H^+ ions on the outer layer of the membrane. As a result, an electrochemical boundary layer with an electrical potential difference forms at this point. An integrated Ag/AgCl reference system serves as the required reference electrode.

The measured voltage is converted to the corresponding pH value using the Nernst equation.

Measuring system



■ 1 Measuring system

- 1 Transmitter CML18
- 2 Smartphone with Smartblue app (optional)
- 3 Memosens CPL57E

Communication and data processing

Communication with the handheld device



Always connect digital laboratory sensors with Memosens technology to a handheld device with Memosens technology, e.g. CML18.

Digital laboratory sensors can store measuring system data in the sensor, including:

- Manufacturer data
 - Serial number
 - Order code
 - Date of manufacture
- Calibration data
 - Calibration date
 - Number of calibrations
 - Serial number of the handheld device used to perform the last calibration or adjustment
- Application data
 - Temperature application range
 - pH application range
 - Date of initial commissioning

Input

Measured variable

pH value

Temperature

Measuring range

- pH value: 0 to 14 pH (1 to 12 pH application range)
- Temperature: -5 to 100 °C (23 to 212 °F) (0 to 80 °C (32 to 176 °F) application range)

Performance characteristics

Reference system

Ag/AgCl lead, bridging electrolyte: gel KCl, 3M, AgCl-free

Environment

Ambient temperature range

NOTICE

Risk of damage from frost!

Do not use the sensor at temperatures below -15 °C (5 °F).

Storage temperature

0 to 50 $^{\circ}$ C (32 to 122 $^{\circ}$ F), non-freezing

Degree of protection

IP 68 (10 m (33 ft) water column, 25 °C (77 °F), 45 days, 1 M KCl)

Electromagnetic compatibility (EMC)

Interference emission and interference immunity as per EN 61326-1: 2013


Process

Process temperature range

0 to 80 °C (32 to 176 °F)

Mechanical construction

Design, dimensions



₽ 2 Engineering unit: mm (in)

Endress+Hauser 3

Weight	40 g (1.4 oz)				
Materials	Sensor shaft	Glass			
	Metal lead	Ag/AgCl			
	Nameplate	Ceramic metal oxide			
	Junction	PTFE			
Temperature sensor	NTC 30K				
Plug-in head	Memosens laboratory plug-in head for digital, non-contact data transmission				
Process connections	Pg 13.5				

Accessories

The following are the most important accessories available at the time this documentation was

▶ For accessories not listed here, please contact your Service or Sales Center.

Device-specific accessories

Memosens data cable CYK10

- For digital sensors with Memosens technology
- Product Configurator on the product page: www.endress.com/cyk10



Technical Information TI00118C

Memosens laboratory cable CYK20

- For digital sensors with Memosens technology
- Product Configurator on the product page: www.endress.com/cyk20

Liquiline Mobile CML18

- Multiparameter mobile device for laboratory and field
- Reliable transmitter with display and app connection
- Product Configurator on the product page: www.endress.com/CML18



Operating Instructions BA02002C

Memobase Plus CYZ71D

- PC software to support laboratory calibration
- Visualization and documentation of sensor management
- Sensor calibrations stored in database
- Product Configurator on the product page: www.endress.com/cyz71d



Technical Information TI00502C

High-quality buffer solutions from Endress+Hauser - CPY20

The secondary buffer solutions have been referenced to primary reference material of the PTB (German Federal Physico-technical Institute) or to standard reference material of NIST (National Institute of Standards and Technology) according to DIN 19266 by a laboratory accredited by the DAkkS (German accreditation body) according to DIN 17025.

Product Configurator on the product page: www.endress.com/cpy20

Endress+Hauser



www.addresses.endress.com