Products

Services

Technical Information Rxn-10 Raman spectroscopic probe





Table of Contents

Function and system design 3
Fields of application3
Laser safety interlock3
Rxn-10 probe3
Rxn-10 probe optics4
Installation5


Specifications	6
Probe specifications	6
Probe dimensions	7
MPE: ocular exposure	8
MPE: skin exposure	8


Function and system design

Fields of application	The Rxn-10 Raman spectroscopic probe is designed for sample measurements in a laboratory, process development, or manufacturing (when part of a single-use probe system) environment. The probe head is compatible with a wide range of interchangeable, commercially available optics (immersion and non-contact) to meet the requirements of different applications. Recommended applications include:				
	 Chemical: reaction monitoring, blending, catalyst monitoring, hydrocarbon speciation, process unit optimization 				
	 Polymer: polymerization reaction monitoring, extrusion monitoring, polymer blending 				
	 Pharmaceutical: active pharmaceutical ingredient (API) reaction monitoring, crystallization Biopharmaceutical: cell culture and fermentation monitoring, optimization, control Food and beverage: zonal heterogeneity mapping of meats and fish 				
	Use of the device for any purpose other than that described, poses a threat to the safety of people and of the entire measuring system and invalidates any warranty.				
Laser safety interlock	The Rxn-10 probe, as installed, forms part of the interlock circuit. If the fiber cable is severed, the laser will turn off within milliseconds of the breakage.				
	NOTICE				
	Handle probes and cables with care.				
	Fiber cables should NOT be kinked and should be routed to maintain the minimum bend radius of 152.4 mm (6 in.).				

• Permanent damage may result if cables are not routed appropriately.

Rxn-10 probe





#	Name	Description
1	Fiber cable	 Connects the probe to the Raman Rxn analyzer via one of the following: Fiber channel (FC) cable assembly Electro-optical (EO) fiber cable
2	Laser emission indicator	When there is potential for the laser to be energized, the indicator light is illuminated.
3	Laser beam shutter	Can be closed to prevent laser emission. Position "I" indicates emission potential. Moving the lever passed position "O" indicates emission is shuttered.
4	Thumb screw	Tighten to secure optics onto the probe when there is not a threaded interface.
5	Optics interface	Insert optics or threaded adapter.

Table 1. Rxn-10 probe parts

Rxn-10 probe optics

The probe is compatible with the following optics to meet the requirements of different applications:

	Optics	Applications
Non-contact optics		For use with solids or turbid media. Also well-suited for delicate or corrosive liquids when sample contamination or damage to optical components is a concern.
Immersion optics (IO)	A0048411	For use in reaction vessels, laboratory reactors, or process streams.
bIO-Optic	A008412	For use with continuous inline measurement in benchtop bioreactor/ fermenter applications requiring headplate entry.
Bio multi optic and bio sleeve	A0051184	For use with continuous inline measurement in benchtop bioreactor/ fermentor applications requiring headplate entry.
Raman optic system for single use	A0048413	For use with disposable fittings for single-use applications.




Table 2. Optics and applications

Installation

During installation, standard eye and skin safety precautions for Class 3B laser products (as per EN-60825/IEC 60825-14 or ANSI Z136.1) should be observed as described below.

	 Standard precautions for laser products should be observed. Probes should always be shuttered or pointed away from people toward a diffuse target if not installed in a sample chamber. 			
	The laser input into the Rxn-10 probe must not exceed 499 mW. If stray light is allowed to enter an unused probe, it will interfere with data collected from a used probe and may cause calibration failure or measurement errors.			
	 Unused probes should ALWAYS be shuttered to prevent stray light from entering the probe. If an optic cap is available, place it on the unused optic. 			
NOTICE	When installing the probe <i>in situ,</i> the user must provide the strain relief to the fiber optic cable at the probe installation location.			

Specifications

Probe specifications

Specifications for the Rxn-10 probe are listed below.

Item		Description		
Laser wavelength	with non-contact or immersion optic	532 nm, 785 nm, or 993 nm		
	with bIO-Optic or Raman optic system for single use	785 nm or 993 nm		
	with bio multi optic and bio sleeve or micro flow bench and micro flow cell	785 nm		
Maximum laser pow	er into probe head	< 499 mW		
Working distance		Based on the sampling optic selected		
Sample interface		Based on the sampling optic selected		
Polarization at samp	le	Unpolarized		
Probe temperature		–10 to 70 °C (14 to 158 °F)		
Temperature ramp		≤ 30 °C/min (≤ 54 °F/min)		
Probe relative humic	lity	20 to 60 %, non-condensing		
Probe spectral cover	age	Probe spectral coverage is limited by the coverage of the analyzer being used		
Laser power at sample	532 nm (with standard 120-mW laser)	> 45 mW		
	785 nm (with standard 400-mW laser)	> 150 mW		
	993 nm (with standard 400-mW laser)	> 150 mW		
Materials of construction	probe body	6061 aluminum, 316L stainless steel, and 303 stainless steel		
	fiber optic cable	Design: PVC jacketed, proprietary construction Connections: proprietary electro-optic or FC to EO fiber converter(s) for non-embedded systems		
Probe	length (not including fiber cable bend radius)	203 mm (8 in.)		
	length (including fiber cable bend radius)	356 mm (14.02 in.)		
	diameter (not including cable)	19 mm (0.75 in.)		
	weight (including cable)	0.5 kg (approximately 1 lb.)		
Fiber optic cable	temperature*	-40 to 70 °C (-40 to 158 °F)		
	length	5 to 25 m (16.4 to 82.0 ft.) lengths standard in 5 m (16.4 ft.) increments		
		Extension fiber cables are also available in lengths from 5 to 200 m (16.4 to 656.2 ft.) in 5 m (16.4 ft.) increments, limited by application.		
	minimum bend radius	152.4 mm (6 in.)		
	flame resistance	Certified: CSA-C/US AWM I/II, A/B, 80C, 30V, FT1, FT2, VW-1, FT4		

* While the fiber optic cable can withstand temperatures up to 80 °C (17 °F), the interface of the cable to the probe head is limited to 70 °C (158 °F).

Table 3. Rxn-10 probe specifications

Probe dimensions

The dimensions for the Rxn-10 probe are shown below.




Figure 2. Rxn-10 probe dimensions

Dimension	Measurement	Description
L1	111 mm 4.37 in.	Length of probe body without cable or optics
L2	203 mm 8 in.	Length with fiber optic cable connected Note: This does not include additional 6 in. minimum bend radius of cable
Н	33 mm 1.3 in.	Height of probe including thumb screw
Ø	19 mm 0.75 in.	Diameter of probe, not including cable

Table 4. Rxn-10 probe dimensions

MPE: ocular exposure

Refer to the tables below from the ANSI Z136.1 standard to calculate the maximum permissible exposure (MPE) for point source ocular exposure to a laser beam.

A correction factor (C_A) may also be required and can be determined below.

Wavelength λ (nm)	Correction factor C _A
400 to 700	1
700 to 1050	10 ^{0.002} (λ-700)
1050 to 1400	5

Table 5. Wavelength dependent correction factor CA

Maximum permissible exposure (MPE) for point source ocular exposure to a laser beam				
Wavelength	Exposure duration t (s)	MPE calculation		
λ (nm)		(J ∙cm⁻²)	(W·cm⁻²)	
532	10 ⁻¹³ to 10 ⁻¹¹	1.0×10^{-7}	-	
	$10^{\text{-}11}$ to 5 \times $10^{\text{-}6}$	2.0×10^{-7}	-	
	5 × 10 ⁻⁶ to 10	1.8 $t^{0.75} \times 10^{-3}$	-	
	10 to 30,000	-	1 × 10 ⁻³	

Table 6. MPE for ocular exposure with 532 nm laser emission

Maximum permissible exposure (MPE) for point source ocular exposure to a laser beam				
Wavelength	Exposure duration	MPE calculation		MPE where
λ (nm)	t (s)	(J·cm⁻²)	(W·cm⁻²)	<i>C</i> _A = 1.4791
785 and 993	10 ⁻¹³ to 10 ⁻¹¹	$1.5 C_{\rm A} \times 10^{-8}$	-	2.2 × 10 ⁻⁸ (J·cm ⁻²)
	10 ⁻¹¹ to 10 ⁻⁹	2.7 $C_{\rm A} t^{0.75}$	-	Insert time (<i>t</i>) and calculate
	10 ⁻⁹ to 18 × 10 ⁻⁶	5.0 <i>C</i> _A × 10 ⁻⁷	-	7.40 × 10 ⁻⁷ (J⋅cm ⁻²)
	18 × 10 ⁻⁶ to 10	$1.8 C_{\rm A} t^{0.75} \times 10^{-3}$	-	Insert time (<i>t</i>) and calculate
	$10 \text{ to } 3 \times 10^4$	-	$C_{\rm A} \times 10^{-3}$	1.4971 × 10 ⁻³ (W·cm ⁻²)

Table 7. MPE for ocular exposure with 785 nm or 993 nm laser emission

MPE: skin exposure

Refer to the table below from the ANSI Z136.1 standard to calculate the MPE for skin exposure to a laser beam.

Maximum permissible exposure (MPE) for skin exposure to a laser beam				
Wavelength	Exposure duration	Exposure duration MPE calcula		MPE where
λ (nm)	t (s)	(J·cm⁻²)	(W·cm ⁻²)	<i>C</i> _A = 1.4791
532, 785 and 993	10 ⁻⁹ to 10 ⁻⁷	2 $C_{\rm A} \times 10^{-2}$	-	2.9582 × 10 ⁻² (J·cm ⁻²)
	10 ⁻⁷ to 10	$1.1 C_{\rm A} t^{0.25}$	-	Insert time (<i>t</i>) and calculate
	10 to 3 × 10 ⁴	-	0.2 <i>C</i> _A	2.9582 × 10 ⁻¹ (W⋅cm ⁻²)

Table 8. MPE for skin exposure with 532 nm, 785 nm or 993 nm laser emission

www.addresses.endress.com

