Endress+Hauser [국]

Technical Information **Memosens CLS82E**

Hygienic conductivity sensor Digital with Memosens technology

Cell constant $k = 0.57 \text{ cm}^{-1}$

Application

For measurements where very diverse conductivities must be measured in one measuring system.

Typical applications include:

- Phase separations
- Chromatography
- Fermentations
- CIP monitoring in small pipes
- Ultrafiltration
- Cleaning of ballast water on ships
- Cleaning of water in a ship's wake

Sensors with temperature probes are used in conjunction with conductivity measuring devices that support automatic temperature compensation:

- Liquiline CM442/CM444/CM448
- Liquiline CM42
- Liquiline CM14

The resistivity in $M\Omega \cdot cm$ can also be measured using these transmitters.

Your benefits

- High measuring accuracy as cell constant is individually measured
- Manufacturer inspection certificate stating the individual cell constant
- Hygienic process connections for installation in pipes or flow vessel
- Easy to clean thanks to electropolished surfaces
- Can be sterilized up to 140 °C (284 °F)
- Stainless steel 1.4435 (AISI 316L) meets the highest demands of the pharmaceutical industry
- IP68 protection
- The entire sensor is certified according to EHEDG and 3-A
- FDA conformity



[Continued from front page]

Other advantages provided by Memosens technology

- Maximum process safetyData security thanks to digital data transmission
- Very easy to use as sensor data are saved in the sensor
- Predictive maintenance can be performed by recording sensor load data in the sensor

Table of contents


Function and system design Measuring principle Measuring system	4 4 4
Communication and data processing	5
Dependability	5 5 5 5
Input Measured variables Measuring ranges Cell constant Temperature compensation	6 6 6 6
Power supply	6 6
Performance characteristics	6
Installation	7 7
Environment	8 9 9 9
Process	9 9 9 9
Mechanical construction	10 10 12 12 12
Certificates and approvals	12
Ordering information Product page Product Configurator Scope of delivery	12 12 12 12

Accessories	12
Measuring cable	13
Calibration solutions	13

Function and system design

Measuring principle


The measuring cell has four electrodes. An alternating current is applied via the outer electrode pair. At the same time, the voltage applied is measured at the two inner electrodes. The electrolytic conductivity between the electrodes can be reliably established based on the measured voltage and the current flow caused by the liquid's resistance. The advantage of this technology compared to traditional two-electrode sensors is that electrochemical effects at the live electrodes are suppressed by the two additional voltage measuring electrodes.



- 1 Conductivity measurement
- I Current intensity measurement
- U Voltage measurement
- G Generator

Measuring system

- A complete measuring system comprises at least:
- Conductivity sensor Memosens CLS82E
- Transmitter, e.g. Liquiline M CM42
- Measuring cable, e.g. Memosens data cable CYK10



- ☑ 2 Example of a measuring system
- 1 Memosens CLS82E
- 2 Transmitter Liquiline CM44x
- 3 Measuring cable

Communication and data processing

Communication with the transmitter

Always connect digital sensors with Memosens technology to a transmitter with Memosens technology. Data transmission to a transmitter for analog sensors is not possible.

Digital sensors can store measuring system data in the sensor. These include the following:


- Manufacturer data
 - Serial number
 - Order code
 - Date of manufacture
- Calibration data
 - Calibration date
 - Cell constant
 - Delta cell constantNumber of calibrations
 - Serial number of the transmitter used to perform the last calibration or adjustment
- Application data
 - Temperature application range
 - Conductivity application range
 - Date of initial commissioning
 - Maximum temperature value
 - Hours of operation at high temperatures

Dependability

Reliability	Memosens technology digitizes the measured values in the sensor and transmits the data to the transmitter via a . The result: If the sensor fails or there is an interruption in the connection between the sensor and transmitter,
	this is reliably detected and reported.The availability of the measuring point is reliably detected and reported.
Maintainability	 Easy handling Sensors with Memosens technology have integrated electronics that store calibration data and other information (e.g. total hours of operation or operating hours under extreme measuring conditions). Once the sensor has been connected, the sensor data are transferred automatically to the transmitter and used to calculate the current measured value. As the calibration data are stored in the sensor, the sensor can be calibrated and adjusted independently of the measuring point. The result: Easy calibration in the measuring lab under optimum external conditions increases the quality of the calibration. Pre-calibrated sensors can be replaced quickly and easily, resulting in a dramatic increase in the availability of the measuring point. Thanks to the availability of the sensor data, maintenance intervals can be accurately defined and predictive maintenance is possible. The sensor history can be documented with external data carriers and evaluation programs. Thus, the current application of the sensors can be made to depend on their previous history.
Integrity	 With inductive transmission of the measured value using a non-contact connection, Memosens guarantees maximum process safety and offers the following benefits: All problems caused by moisture are eliminated. Plug-in connection remains free from corrosion Measured value distortion from moisture is not possible. The plug-in system can even be connected under water. The transmitter is galvanically decoupled from the medium. EMC safety is guaranteed by screening measures for the digital transmission of measured values.

	Input	
Measured variables	ConductivityTemperature	
Measuring ranges	Conductivity ¹⁾	1 µS/cm to 500 mS/cm
	1) In relation to water at 25 $^{\circ}$ C (77 $^{\circ}$ F)	
	Temperature	-5 to 140 °C (23 to 284 °F)
Cell constant	k = 0.57 cm ⁻¹	
Temperature compensation	Pt1000 (Class A according to IEC 6075	1)

Power supply



Performance characteristics

Measuring uncertainty	measuring system traceable to NIST of	sured in a solution with approx. 50 μS/cm using a reference or PTB. The exact cell constant is entered into the manufacturer ncertainty of measurement in determining the cell constant is	
Response time	Conductivity	t ₉₅ ≤ 2 s	
	Temperature ¹⁾		
	With Pg 13.5 or Clamp	$t_{90} \le 16 \text{ s}^{2)}$	
	With other process connection	$t_{90} \le 28 \ s^{2)}$	
	 DIN VDI/VDE 3522-2 (0.3 m/s laminar) With temperature prediction activated as standard 		
Measured error	Conductivity		
	In the range 1 μ S/cm to 1 mS/cm 1	$^{)} \leq 2$ % of reading	
	In the range 1 mS/cm to 500 mS/c	$m^{(1)} \le 4\%$ of reading	
	Temperature		
	With Pg 13.5 or Clamp	≤ 0.5 K, in measuring range -5 to 100 °C (23 to 212 °F) ≤ 1.0 K, in measuring range 100 to 140 °C (212 to 284 °F)	
	With other process connection	\leq 1.0 K, in measuring range -5 to 140 °C (23 to 284 °F)	
	1) In as-delivered state (factory adjust	ment at 50 μS/cm)	


Repeatability

Conductivity Temperature \leq 0.2 % of reading, in specified measuring range \leq 0.05 K

Installation

Installation instructions

Symmetrical installation is recommended in order to guarantee linearity. The distance to the side walls and opposite walls must be at least 15 mm.




Minimum distance between pipe and end of the measuring cell

- 1 Sensor
- 2 Pipe
- 3 Direction of flow


The ionic current in the liquid is affected by the walls in confined installation conditions. This effect is compensated by what is referred to as the installation factor. The installation factor can be entered in the transmitter for the measurement or the cell constant is corrected by multiplying by the installation factor.

The value of the installation factor depends on the diameter and the conductivity of the pipe nozzle as well as the sensor's distance to the wall. The installation factor can be disregarded (f = 1.00) if the distance to the wall is sufficient (a > 15 mm). If the distance to the wall is smaller, the installation factor increases for electrically insulating pipes (f > 1) and decreases for electrically conductive pipes (f < 1). The installation factor can be determined using calibration solutions.



I 5 Schematic drawing of the sensor in confined installation conditions


- a Wall distance
- b Gap width



Relationship between installation factor f and wall distance a

1 Electrically insulating pipe wall

2 Electrically conductive pipe wall



Relationship between installation factor f and gap width b

- 1 Electrically insulating pipe wall
- 2 Electrically conductive pipe wall

Hygienic properties

For 3-A-compliant installation, please observe the following:


- After the device has been mounted, hygienic integrity must be guaranteed.
- ► 3-A-compliant process connections must be used.

Installation factors for assemblies

Environment

Ambient temperature


−20 to 60 °C (−4 to 140 °F)



8 Pressure/temperature ratings

A Can be sterilized for a short time (45 min.)

A0044758




Mechanical construction

Dimensions in mm (in)

🖻 11 DN 25 brown

🖻 12 DN 25 standard



Weight	Depending on the version, e.g. Process connection Pg 13.5: 0.06 to 0.09 kg (0.13 to 0.20 lbs) Process connection G1 or NPT: approx. 0.9 kg (1.98 lbs)	
Materials (in contact with medium)	Sensor element: Platinum and ceramic (zirconium oxide) Process connection: Stainless steel 1.4435 (AISI 316L) Only for CLS82E-**NA* ¹⁾ and CLS82E-**NB* ²⁾ : Seal: EPDM 1) Connection DN25 standard 2) Connection DN25 B. Braun	
Surface roughness	$R_a < 0.38 \ \mu m$	

Certificates and approvals

Current certificates and approvals for the product are available via the Product Configurator at www.endress.com.

1. Select the product using the filters and search field.

2. Open the product page.

The **Configuration** button opens the Product Configurator.

Ordering information

Product page	www.endress.com/cls82e
Product Configurator	 On the product page there is a Configure button to the right of the product image. 1. Click this button. The Configurator opens in a separate window.
	 Select all the options to configure the device in line with your requirements. In this way, you receive a valid and complete order code for the device. Export the order code as a PDF or Excel file. To do so, click the appropriate button on the right above the selection window.
	For many products you also have the option of downloading CAD or 2D drawings of the selected product version. Click the CAD tab for this and select the desired file type using picklists.
Scope of delivery	The scope of delivery includes: Sensor in the version ordered Operating Instructions

Accessories

The following are the most important accessories available at the time this documentation was issued.

► For accessories not listed here, please contact your Service or Sales Center.

Measuring cable	 Memosens data cable CYK10 For digital sensors with Memosens technology Product Configurator on the product page: www.endress.com/cyk10
	Technical Information TI00118C
	 Memosens data cable CYK11 Extension cable for digital sensors with Memosens protocol Product Configurator on the product page: www.endress.com/cyk11 Technical Information TI00118C
Calibration solutions	Conductivity calibration solutions CLY11
	 Precision solutions referenced to SRM (Standard Reference Material) by NIST for qualified calibration of conductivity measuring systems in accordance with ISO 9000 CLY11-A, 74 µS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081902 CLY11-B, 149.6 µS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081903 CLY11-C, 1.406 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081904 CLY11-D, 12.64 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081905 CLY11-E, 107.00 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz) Order No. 50081905
	Technical Information TI00162C



www.addresses.endress.com

