Technical Information Fiber Optic Cables





Technical Information Fiber Optic Cables

Table of Contents

Function and system design		
Introduction	3	
Raman cables	3	

Specifications	4
General	4
Cable types	4

Fiber Optic Cables Technical Information

Function and system design

Introduction

Fiber optic technology revolutionized Raman spectroscopy by allowing Raman sampling probes to be located remote from a base unit. This enabled Raman spectra to be acquired in hazardous environments from samples that cannot be easily transported to a sampling chamber. Consequently, Raman spectroscopy entered several new arenas including the industrial process line, where the base unit is placed in a control room or other protected environment while the Raman probe is placed in the process line for real-time, *in situ* process monitoring and control.

In the majority of state-of-the-art remote dispersive Raman systems, the excitation radiation is delivered from the laser to the Raman probe through a single excitation fiber. The scattered radiation that is collected from the sample is delivered to the spectrograph through a single collection fiber.

Optical fibers are constructed of a low-hydroxyl silica core surrounded by a flourine-doped silica cladding and a protective acrylate buffer coating (this 3-layer fiber is typically formed in a single "draw" manufacturing operation). The outer packaging of the cable may vary depending on the application. Fibers intended for industrial and laboratory applications often place a tight polymer buffer on the fiber and/or run through a loose polymer tube. Such fiber subassemblies can then be packaged into a composite industrial grade cable with a robust polymer outer jacket containing other such optical fiber subassemblies, electrical wires, and a rigid strength member.

Raman cables

All Endress+Hauser Raman probes use standard cables comprised of an integrated fiber cable assembly containing an excitation fiber and a collection fiber packaged in a robust Polyvinyl Chloride (PVC) jacket to prevent breakage. Endress+Hauser Raman fiber optic probes also have integrated the laser interlock into the probe termination for improved laser safety. The laser switches off within milliseconds, preventing laser light from being dispersed into the environment if the cable is severed.

Endress+Hauser Raman's standard fiber cables are indoor/outdoor, riser-rated fiber optic cables. These cables are also fully rated for flame/UV resistance and pull strength, maximizing their safety in the process environment. Endress+Hauser Raman fiber cables are suitable for use in a variety of environments including direct burial, underground ducts, aerial installations, steam tunnels, building risers, cable trays, and harsh industrial settings.




Figure 1. Cross section of a Raman fiber optic cable

Technical Information Fiber Optic Cables

Specifications

General

The specifications for fiber optic cables are listed below:

Item	Description
General Features	Integrated copper conductor wire for interlock capability Kevlar internal strength members Flame retardant Fungus resistant
Cable Rating	Operating temperature: -40°C to +80°C Storage temperature: -55°C to +80°C Indoor/outdoor UV Riser-rated Cable tray-rated Certified: CSA-C/US AWM I/II, A/B, 80C, 30V, FTI, FT2, VW-1, FT4 Rated: AWM I/II A/B 80C 30V FT4
Bend Radius	Installation: 6.3" (16.00 cm) Operation: 3.2" (8.13 cm)
Crush Resistance	1700-2200 N/cm
Termination	Proprietary electro-optic

Table 1. Specifications

Cable types

Fiber optic cables with different connectors are available to connect various Raman probes and Raman Rxn analyzers. A list of commonly used fiber cables is provided below.

2011654-XXX



Figure 2. 2011654-XXX

Analyzer	Probe	Description	Standard Length
Raman Rxn5	Rxn- 30	FOCA, EO(M)/SSCS, XXXMSpecification: Electro/Optical Cable; Length = XXX Meters; Connectors = EO (M) to SSCS	No standard length (limited by application)

Table 2. 2011654-XXX

2018150-XXX



Figure 3. 2018150-XXX

Analyzer	Probe	Description	Standard Length
Raman Rxn2, Raman Rxn4, Raman Rxn5, Legacy Rxn products	Probes that accept EO connectors; Probes that accept FC connectors	FOCA, EO(M)/FC, XXXM Specification: Electro/Optical Cable; Length = XXX Meters; Connectors = EO (M) to FC	5-200 meters in 5-meter increments (limited by application)

Table 3. 2018150-XXX

NOTICE

► This fiber optic cable is compatible with some legacy Rxn products.

Fiber Optic Cables Technical Information

2017161-XXX



Figure 4. 2017161-XXX

Analyzer	Probe	Description	Standard Length
Raman Rxn2, Raman Rxn4, Raman Rxn5	Probes that accept EO connectors	FOCA, EO(M)/EO(M), XXXM Specification: Electro/Optical Cable; Length = XXX Meters; Connectors = EO (M) to EO (M)	5-200 meters in 5-meter increments (limited by application)

Table 4. 2017161-XXX (used as an extension cable for 2012936)

2012936-XXX



Figure 5. 2012936-XXX

Analyzer	Probe	Description	Standard Length
Raman Rxn2, Raman Rxn4, Raman Rxn5	Probes that accept EO connectors	FOCA, EO(M)/EO(M), XXXM Specification: Electro/Optical Cable; Length = XXX Meters; Connectors = EO (M) to EO (M)	5-200 meters in 5-meter increments (limited by application)

Table 5. 2012936-XXX

2018539-XXX



Figure 6. 2018539-XXX

Analyzer	Probe	Description	Standard Length
Legacy Rxn platforms	Probes that accept FC connectors	FOCA, EXT, FC/FC, XXXMSpecification: Electro/Optical Cable; Length = XXX Meters; Connectors = FC to FC	5-200 meters in 5-meter increments (limited by application)

Table 6. 2018539-XXX

NOTICE

▶ This fiber optic cable is compatible with some legacy Rxn products.

2018540-XXX



Figure 7. 2018540-XXX

Analyzer	Probe	Description	Standard Length
Legacy Rxn platforms	Probes that accept FC connectors	Electro/Optical Cable; Length = XXX	5-200 meters in 5-meter increments (limited by application)

Table 7. 2018540-XXX

NOTICE

This fiber optic cable is compatible with some legacy Rxn products.

Technical Information Fiber Optic Cables

www.addresses.endress.com



People for Process Automation