

# Operating Instructions Indumax CLS54D

Hygienic inductive conductivity sensor  
For the food and beverage industry, the  
pharmaceutical industry and biotechnology









# Table of contents








<b>1</b>	<b>About this document</b> .....	<b>4</b>	<b>Index</b> .....	<b>20</b>
1.1	Warnings .....	4		
1.2	Symbols .....	4		
1.3	Symbols on the device .....	4		
1.4	Documentation .....	5		
<b>2</b>	<b>Basic safety instructions</b> .....	<b>5</b>		
2.1	Requirements for the personnel .....	5		
2.2	Intended use .....	5		
2.3	Workplace safety .....	6		
2.4	Operational safety .....	6		
2.5	Product safety .....	6		
<b>3</b>	<b>Incoming acceptance and product identification</b> .....	<b>6</b>		
3.1	Incoming acceptance .....	6		
3.2	Product identification .....	7		
3.3	Scope of delivery .....	8		
<b>4</b>	<b>Mounting</b> .....	<b>8</b>		
4.1	Mounting requirements .....	8		
4.2	Mounting the sensor .....	11		
4.3	Post-mounting check .....	11		
<b>5</b>	<b>Electrical connection</b> .....	<b>12</b>		
5.1	Connecting the sensor .....	12		
5.2	Ensuring the degree of protection .....	13		
5.3	Post-connection check .....	13		
<b>6</b>	<b>Maintenance</b> .....	<b>14</b>		
<b>7</b>	<b>Repair</b> .....	<b>15</b>		
7.1	General notes .....	15		
7.2	Spare parts .....	15		
7.3	Return .....	15		
7.4	Disposal .....	15		
<b>8</b>	<b>Accessories</b> .....	<b>16</b>		
8.1	Cable extension .....	16		
8.2	Calibration solutions .....	16		
<b>9</b>	<b>Technical data</b> .....	<b>17</b>		
9.1	Input .....	17		
9.2	Performance characteristics .....	17		
9.3	Environment .....	17		
9.4	Process .....	18		
9.5	Mechanical construction .....	19		

# 1 About this document



## 1.1 Warnings

Structure of information	Meaning
<p> <b>DANGER</b></p> <p><b>Causes (/consequences)</b> If necessary, Consequences of non-compliance (if applicable)</p> <ul style="list-style-type: none"> <li>▶ Corrective action</li> </ul>	<p>This symbol alerts you to a dangerous situation. Failure to avoid the dangerous situation <b>will</b> result in a fatal or serious injury.</p>
<p> <b>WARNING</b></p> <p><b>Causes (/consequences)</b> If necessary, Consequences of non-compliance (if applicable)</p> <ul style="list-style-type: none"> <li>▶ Corrective action</li> </ul>	<p>This symbol alerts you to a dangerous situation. Failure to avoid the dangerous situation <b>can</b> result in a fatal or serious injury.</p>
<p> <b>CAUTION</b></p> <p><b>Causes (/consequences)</b> If necessary, Consequences of non-compliance (if applicable)</p> <ul style="list-style-type: none"> <li>▶ Corrective action</li> </ul>	<p>This symbol alerts you to a dangerous situation. Failure to avoid this situation can result in minor or more serious injuries.</p>
<p> <b>NOTICE</b></p> <p><b>Cause/situation</b> If necessary, Consequences of non-compliance (if applicable)</p> <ul style="list-style-type: none"> <li>▶ Action/note</li> </ul>	<p>This symbol alerts you to situations which may result in damage to property.</p>

## 1.2 Symbols

	Additional information, tips
	Permitted or recommended
	Not permitted or not recommended
	Reference to device documentation
	Reference to page
	Reference to graphic
	Result of a step

## 1.3 Symbols on the device

	Reference to device documentation
	Do not dispose of products bearing this marking as unsorted municipal waste. Instead, return them to the manufacturer for disposal under the applicable conditions.

## 1.4 Documentation

The following manuals, which complement these Operating Instructions, can be found on the product pages on the Internet:



Technical Information Indumax CLS54D, TI00508C

# 2 Basic safety instructions

## 2.1 Requirements for the personnel

- Installation, commissioning, operation and maintenance of the measuring system may be carried out only by specially trained technical personnel.
- The technical personnel must be authorized by the plant operator to carry out the specified activities.
- The electrical connection may be performed only by an electrical technician.
- The technical personnel must have read and understood these Operating Instructions and must follow the instructions contained therein.
- Faults at the measuring point may only be rectified by authorized and specially trained personnel.



Repairs not described in the Operating Instructions provided must be carried out only directly at the manufacturer's site or by the service organization.

## 2.2 Intended use

Indumax CLS54D is designed for the inductive measurement of the conductivity of liquids in the food and beverages industry.

The six-decade measuring range and the excellent chemical resistance properties of the materials in contact with the medium make it possible to use this sensor in a wide range of applications, such as:

- Concentration measurement of acids and bases
- Phase separation of products

The sensor is used with Liquiline CM44x/R/P or Liquiline CM42.

Use of the device for any purpose other than that described, poses a threat to the safety of people and of the entire measuring system and is therefore not permitted.

The manufacturer is not liable for damage caused by improper or non-designated use.

### NOTICE

#### Applications outside specifications!

Incorrect measurements, malfunctions and even measuring point failure could result

- ▶ Use the product only in accordance with the specifications.
- ▶ Pay attention to the technical data on the nameplate.

## 2.3 Workplace safety

As the user, you are responsible for complying with the following safety conditions:

- Installation guidelines
- Local standards and regulations

### Electromagnetic compatibility

- The product has been tested for electromagnetic compatibility in accordance with the applicable international standards for industrial applications.
- The electromagnetic compatibility indicated applies only to a product that has been connected in accordance with these Operating Instructions.

## 2.4 Operational safety

### Before commissioning the entire measuring point:

1. Verify that all connections are correct.
2. Ensure that electrical cables and hose connections are undamaged.
3. Do not operate damaged products, and protect them against unintentional operation.
4. Label damaged products as defective.

### During operation:

- ▶ If faults cannot be rectified:  
products must be taken out of service and protected against unintentional operation.

## 2.5 Product safety

The product is designed to meet state-of-the-art safety requirements, has been tested, and left the factory in a condition in which it is safe to operate. The relevant regulations and international standards have been observed.

# 3 Incoming acceptance and product identification

## 3.1 Incoming acceptance

1. Verify that the packaging is undamaged.
  - ↳ Notify the supplier of any damage to the packaging.  
Keep the damaged packaging until the issue has been resolved.
2. Verify that the contents are undamaged.
  - ↳ Notify the supplier of any damage to the delivery contents.  
Keep the damaged goods until the issue has been resolved.
3. Check that the delivery is complete and nothing is missing.
  - ↳ Compare the shipping documents with your order.

4. Pack the product for storage and transportation in such a way that it is protected against impact and moisture.
  - ↳ The original packaging offers the best protection.  
Make sure to comply with the permitted ambient conditions.

If you have any questions, please contact your supplier or your local Sales Center.

## 3.2 Product identification

### 3.2.1 Nameplate

The nameplate provides you with the following information on your device:

- Manufacturer identification
  - Extended order code
  - Serial number
- ▶ Compare the information on the nameplate with the order.

### 3.2.2 Product identification

#### Product page

[www.endress.com/cls54D](http://www.endress.com/cls54D)

#### Interpreting the order code

The order code and serial number of your product can be found in the following locations:

- On the nameplate
- In the delivery papers

#### Obtaining information on the product

1. Go to [www.endress.com](http://www.endress.com).
2. Page search (magnifying glass symbol): Enter valid serial number.
3. Search (magnifying glass).
  - ↳ The product structure is displayed in a popup window.
4. Click the product overview.
  - ↳ A new window opens. Here you fill information pertaining to your device, including the product documentation.

#### Manufacturer address

Endress+Hauser Conducta GmbH+Co. KG  
Dieselstraße 24  
D-70839 Gerlingen

### 3.3 Scope of delivery

The scope of delivery includes:

- Sensor in the version ordered
- Operating Instructions
- ▶ If you have any queries:  
Please contact your supplier or local sales center.

## 4 Mounting

### 4.1 Mounting requirements

#### 4.1.1 Hygienic requirements

- ▶ Easily cleanable installation of equipment according to the criteria of the EHEDG must be free of dead legs.
- ▶ If a dead leg is unavoidable, it shall be kept as short as possible. Under no circumstances shall the length of a dead leg  $L$  exceed the pipe's inner diameter  $D$  minus the equipment's enveloping diameter  $d$ . The condition  $L \leq D - d$  applies.
- ▶ Furthermore, the dead leg must be self-draining, so neither product nor process fluids are retained therein.
- ▶ Within tank installations, the cleaning device must be located so that it directly flushes the dead leg.
- ▶ For further reference, see the recommendations concerning hygienic seals and installations in EHEDG Doc. 10 and the Position Paper: "Easy cleanable Pipe couplings and Process connections".


For 3-A-compliant installation, please observe the following:

- ▶ After the device has been mounted, hygienic integrity must be guaranteed.
- ▶ 3-A-compliant process connections must be used.

#### 4.1.2 Orientation

The sensor must be completely immersed in the medium. Avoid air bubbles in the area of the sensor.





A0037970

**1** Installation positions of the conductivity sensor

**i** Changes in the flow direction (after pipe bends) can cause turbulence in the medium. Install the sensor at a distance of at least 1 m (3.3 ft) downstream from a pipe bend.

The product should flow along the hole of the sensor (see the arrows on the housing). The symmetrical measuring channel allows flow in both directions.

### 4.1.3 Installation factor


The ionic current in the liquid is affected by the walls in confined installation conditions. This effect is compensated by what is referred to as the installation factor. The installation factor can be entered in the transmitter for the measurement or the cell constant is corrected by multiplying by the installation factor.

The value of the installation factor depends on the diameter and the conductivity of the pipe nozzle as well as the distance  $a$  between the sensor and the wall.

The installation factor  $f$  ( $f = 1.00$ ) can be disregarded if the distance to the wall is sufficient ( $a > 15$  mm, from DN 65).


If the distance to the wall is smaller, the installation factor increases for electrically insulating pipes ( $f > 1$ ), and decreases for electrically conductive pipes ( $f < 1$ ).

It can be measured using calibration solutions, or a close approximation can be determined from the diagram below.



2 Installation CLS54D

$a$  Wall distance




3 Relationship between installation factor  $f$  and wall distance  $a$

- 1 Electrically conductive pipe wall
- 2 Electrically insulating pipe wall

### 4.1.4 Air set

The digital sensor has already been adjusted at the factory. Onsite compensation is not required.

## 4.2 Mounting the sensor



A0032586

4 Installed length of the sensor

1 Direction of medium flow

a Distance from pipe wall

- ▶ When installing, align the sensor in such a way that the medium flows through the flow opening of the sensor in the direction of medium flow.
  - ↳ The sensor head must be completely immersed in the medium.

## 4.3 Post-mounting check

Put the sensor into operation only if you can answer yes to the following questions:

1. Are the sensor and cable undamaged?
2. Is the orientation correct?
3. Has the sensor been installed in the process connection, and does not suspend freely from the cable?

## 5 Electrical connection

### **⚠ WARNING**

#### **Device is live!**


Incorrect connection may result in injury or death!

- ▶ The electrical connection may be performed only by an electrical technician.
- ▶ The electrical technician must have read and understood these Operating Instructions and must follow the instructions contained therein.
- ▶ **Prior** to commencing connection work, ensure that no voltage is present on any cable.

### 5.1 Connecting the sensor

The sensor is supplied with a fixed cable. The wiring diagram is provided in the Operating Instructions of the transmitter used.

Connection via a junction box is necessary for a cable connection. The extension to the transmitter is via the CYK11 cable.




A0032587

■ 5 Connection with CYK11 cable extension via junction box, dimensions in mm (inch)

- 1 Cable glands - shield fixed in gland
- 2 Shielding
- A CYK11 from transmitter
- B Sensor cable

Sensors with a fixed cable and M12 plug can be extended with the CYK11 measuring cable and an M12 socket.



A0017842

**6** *CYK11 for extension with M12 connection*

- 1 Transmitter
- 2 CYK11 measuring cable with M12 connection
- A CLS54D connecting cable with M12 plug
- B Sensor CLS54D

## 5.2 Ensuring the degree of protection

Only the mechanical and electrical connections which are described in these instructions and which are necessary for the required, designated use, may be carried out on the device delivered.

- ▶ Exercise care when carrying out the work.

Otherwise, the individual types of protection (Ingress Protection (IP), electrical safety, EMC interference immunity) agreed for this product can no longer be guaranteed due, for example to covers being left off or cable (ends) that are loose or insufficiently secured.

## 5.3 Post-connection check

Device health and specifications	Action
Is the outside of the sensor, assembly or cable free from damage?	▶ Perform a visual inspection.
Electrical connection	Action
Are the mounted cables strain-relieved and not twisted?	<ul style="list-style-type: none"> <li>▶ Perform a visual inspection.</li> <li>▶ Untwist the cables.</li> </ul>
Is a sufficient length of the cable cores stripped, and are the cores positioned in the terminal correctly?	<ul style="list-style-type: none"> <li>▶ Perform a visual inspection.</li> <li>▶ Pull gently to check they are seated correctly.</li> </ul>
Are all screw terminals tightened?	▶ Tighten the screw terminals.
Are all cable entries mounted, firmly tightened and leak-tight?	▶ Perform a visual inspection.
Are all cable entries mounted on the side or pointing downwards?	In the case of lateral cable entries: ▶ Point cable loops downward so that water can drip off.

## 6 Maintenance

### WARNING

#### Thiocarbamide

Harmful if swallowed! Limited evidence of carcinogenicity! Possible risk of harm to the unborn child! Dangerous for the environment with long-term effects!

- ▶ Wear protective goggles, protective gloves and appropriate protective clothing.
- ▶ Avoid all contact with the eyes, mouth and skin.
- ▶ Avoid discharge into the environment.

### CAUTION

#### Corrosive chemicals

Risk of chemical burns to the eyes and skin and risk of damage to clothing and equipment!

- ▶ It is absolutely essential to protect the eyes and hands properly when working with acids, alkalis and organic solvents!
- ▶ Wear protective goggles and safety gloves.
- ▶ Clean away splashes on clothes and other objects to prevent any damage.
- ▶ Comply with instructions in the safety data sheets for the chemicals used.

Clean away fouling on the sensor as follows depending on the type of fouling:

#### 1. Oily and greasy films:

Clean with a grease remover, e.g. alcohol, or hot water and a surfactant-containing (basic) agent (e.g. washing-up liquid).

#### 2. Lime and metal hydroxide buildup and low solubility (lyophobic) organic buildup:

Dissolve buildup with diluted hydrochloric acid (3 %) and then rinse thoroughly with plenty of clear water.

#### 3. Sulfidic buildup (from flue gas desulfurization or wastewater treatment plants):

Use a mixture of hydrochloric acid (3 %) and thiocarbamide (commercially available) and then rinse thoroughly with plenty of clear water.

#### 4. Buildup containing protein (e.g. in the food industry):

Use a mixture of hydrochloric acid (0.5 %) and pepsin (commercially available) and then rinse thoroughly with plenty of clear water.

#### 5. Readily soluble biological buildup:

Rinse with pressurized water.

After cleaning, rinse the sensor thoroughly with plenty of water.

## 7 Repair

### 7.1 General notes

The repair and conversion concept provides for the following:

- The product has a modular design
- Spare parts are grouped into kits which include the associated kit instructions
- Only use original spare parts from the manufacturer
- Repairs are carried out by the manufacturer's Service Department or by trained users
- Certified devices can only be converted to other certified device versions by the manufacturer's Service Department or at the factory
- Observe applicable standards, national regulations, Ex documentation (XA) and certificates

1. Carry out the repair according to the kit instructions.
2. Document the repair and conversion and enter, or have entered, in the Life Cycle Management tool (W@M).

### 7.2 Spare parts

Device spare parts that are currently available for delivery can be found on the website:

[www.endress.com/device-viewer](http://www.endress.com/device-viewer)

- ▶ Quote the serial number of the device when ordering spare parts.

### 7.3 Return

The product must be returned if repairs or a factory calibration are required, or if the wrong product was ordered or delivered. As an ISO-certified company and also due to legal regulations, Endress+Hauser is obliged to follow certain procedures when handling any returned products that have been in contact with medium.

To ensure the swift, safe and professional return of the device:

- ▶ Refer to the website [www.endress.com/support/return-material](http://www.endress.com/support/return-material) for information on the procedure and conditions for returning devices.

### 7.4 Disposal



If required by the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE), the product is marked with the depicted symbol in order to minimize the disposal of WEEE as unsorted municipal waste. Do not dispose of products bearing this marking as unsorted municipal waste. Instead, return them to the manufacturer for disposal under the applicable conditions.

## 8 Accessories

The following are the most important accessories available at the time this documentation was issued.

- For accessories not listed here, please contact your Service or Sales Center.

### 8.1 Cable extension

#### 8.1.1 Measuring cable

##### Memosens data cable CYK11

- Extension cable for digital sensors with Memosens protocol
- Product Configurator on the product page: [www.endress.com/cyk11](http://www.endress.com/cyk11)



Technical Information TI00118C

#### 8.1.2 Junction box

##### Junction box, M12 socket/cable

- Material: aluminum, painted
- Cable extension: Memosens sensors, Liquiline
- Order No.: 71145498

##### Junction box, cable/cable

- Material: aluminum, painted
- Cable extension: Memosens sensors, Liquiline
- Order No.: 71145499

## 8.2 Calibration solutions

### Conductivity calibration solutions CLY11

Precision solutions referenced to SRM (Standard Reference Material) by NIST for qualified calibration of conductivity measuring systems in accordance with ISO 9000

- CLY11-B, 149.6  $\mu\text{S}/\text{cm}$  (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)  
Order No. 50081903
- CLY11-C, 1.406 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)  
Order No. 50081904
- CLY11-D, 12.64 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)  
Order No. 50081905
- CLY11-E, 107.00 mS/cm (reference temperature 25 °C (77 °F)), 500 ml (16.9 fl.oz)  
Order No. 50081906



Technical Information TI00162C



## 9 Technical data

### 9.1 Input

#### 9.1.1 Measured variables

- Conductivity
- Temperature

#### 9.1.2 Measuring range

Conductivity	Recommended range: 100 $\mu\text{S}/\text{cm}$ to 2000 $\text{mS}/\text{cm}$ (uncompensated)
Temperature	-10 to +150 $^{\circ}\text{C}$ (+14 to +302 $^{\circ}\text{F}$ )

#### 9.1.3 Cell constant

$$k = 6.3 \text{ cm}^{-1}$$

#### 9.1.4 Temperature measurement

Pt1000 (Class A according to DIN EN 60751)

### 9.2 Performance characteristics

#### 9.2.1 Conductivity response time

$$t_{95} \leq 2 \text{ s}$$

#### 9.2.2 Temperature response time

$$t_{90} \leq 26 \text{ s}$$

#### 9.2.3 Maximum measured error

< 100 $^{\circ}\text{C}$ (212 $^{\circ}\text{F}$ ):	$\pm(10 \mu\text{S}/\text{cm} + 0.5 \% \text{ of reading})$ , after calibration
> 100 $^{\circ}\text{C}$ (212 $^{\circ}\text{F}$ ):	$\pm(25 \mu\text{S}/\text{cm} + 0.5 \% \text{ of reading})$ , after calibration

#### 9.2.4 Repeatability

$$0.2 \% \text{ of reading} + 3 \mu\text{S}/\text{cm}$$

### 9.3 Environment

#### 9.3.1 Ambient temperature

$$-20 \text{ to } 60 \text{ }^{\circ}\text{C} \text{ (-4 to } 140 \text{ }^{\circ}\text{F)}$$

#### 9.3.2 Storage temperature

$$-25 \text{ to } +80 \text{ }^{\circ}\text{C} \text{ (-13 to } +176 \text{ }^{\circ}\text{F)}$$

### 9.3.3 Relative humidity

5 to 95 %

### 9.3.4 Degree of protection

IP 68 / NEMA type 6P (1 m water column, 25 °C, 168 h)

## 9.4 Process

### 9.4.1 Process temperature

-10 to +125 °C (+14 to +257 °F)

### 9.4.2 Sterilization

150 °C (302 °F) / 6 bar (87 psi) absolute, (max. 60 min.)


### 9.4.3 Process pressure (absolute)

13 bar (188.5 psi) up to 90 °C (194 °F)

9 bar (130.5 psi) at 125 °C (257 °F)

Underpressure down to 0.1 bar (1.45 psi)

### 9.4.4 Pressure/temperature ratings



A0008379

#### 7 Pressure/temperature ratings

A Temporarily for sterilization (max. 60 min.)

B MAWP (maximum allowable working pressure) according to ASME-BPVC Sec. VIII, Div 1 UG101 for CRN registration

## 9.5 Mechanical construction

### 9.5.1 Weight

0.3 to 0.5 kg (0.66 to 1.1 lb.) depending on version plus cable

### 9.5.2 Materials

In contact with medium

Virgin PEEK

Not in contact with medium

PPS-GF40

SMS coupling: stainless steel 1.4301 (AISI 304) or 1.4307 (AISI 304L)

Sanitary coupling: stainless steel 1.4404 (AISI 316L)

Cable gland: PEEK

Seals: FKM,

Cable: TPE

### 9.5.3 Surface roughness

$R_a \leq 0.8 \mu\text{m}$  (smooth, injection-molded PEEK surface) at surfaces in contact with medium

### 9.5.4 Chemical resistance

Medium	Concentration	PEEK
Caustic soda NaOH	0 to 15 %	20 to 90 °C (68 to 194 °F)
Nitric acid HNO <sub>3</sub>	0 to 10 %	20 to 90 °C (68 to 194 °F)
Phosphoric acid H <sub>3</sub> PO <sub>4</sub>	0 to 15 %	20 to 80 °C (68 to 176 °F)
Sulfuric acid H <sub>2</sub> SO <sub>4</sub>	0 to 30 %	20 °C (68 °F)
Peracetic acid H <sub>3</sub> C-CO-OOH	0.2 %	20 °C (68 °F)

# Index

## A

Accessories . . . . .	16
Air set . . . . .	10
Ambient temperature . . . . .	17

## C

Calibration solutions . . . . .	16
Cell constant . . . . .	17
Check	
Connection . . . . .	13
Mounting . . . . .	11
Chemical resistance . . . . .	19
Cleaning agent . . . . .	14
Conductivity response time . . . . .	17
Connection	
Check . . . . .	13
Ensuring the degree of protection . . . . .	13

## D

Degree of protection . . . . .	18
Ensuring . . . . .	13
Disposal . . . . .	15

## E

Electrical connection . . . . .	12
Environment . . . . .	17

## I

Incoming acceptance . . . . .	6
Input . . . . .	17
Installation factor . . . . .	10
Intended use . . . . .	5
Interpreting the order code . . . . .	7

## J

Junction box . . . . .	16
------------------------	----

## M

Maintenance . . . . .	14
Manufacturer address . . . . .	7
Materials . . . . .	19
Maximum measured error . . . . .	17
Measured variables . . . . .	17
Measuring cable . . . . .	16
Measuring ranges . . . . .	17
Mechanical construction . . . . .	19

Mounting . . . . .	8
Mounting requirements . . . . .	8

## N

Nameplate . . . . .	7
---------------------	---

## O

Operational safety . . . . .	6
Orientation . . . . .	8

## P

Performance characteristics . . . . .	17
Post-mounting check . . . . .	11
Pressure/temperature ratings . . . . .	18
Process . . . . .	18
Process pressure . . . . .	18
Process temperature . . . . .	18
Product identification . . . . .	6, 7
Product page . . . . .	7
Product safety . . . . .	6

## R

Relative humidity . . . . .	18
Repair . . . . .	15
Repeatability . . . . .	17
Requirements for the personnel . . . . .	5
Return . . . . .	15

## S

Safety instructions . . . . .	5
Scope of delivery . . . . .	8
Sensor	
Connecting . . . . .	12
Mounting . . . . .	11
Spare parts . . . . .	15
Sterilization . . . . .	18
Storage temperature . . . . .	17
Surface roughness . . . . .	19
Symbols . . . . .	4

## T

Technical data . . . . .	17
Environment . . . . .	17
Mechanical construction . . . . .	19
Performance characteristics . . . . .	17
Process . . . . .	18

---

Temperature measurement . . . . .	17
Temperature response time . . . . .	17
Temperature/pressure ratings . . . . .	18

**U**

Use . . . . .	5
---------------	---

**W**

Warnings . . . . .	4
Weight . . . . .	19
Wiring . . . . .	12
Workplace safety . . . . .	6







71559947

[www.addresses.endress.com](http://www.addresses.endress.com)

---