


Installation Instructions

Adapter RSE8, RSE5, RSE4

Proline Dosimass, Dosimag

EN Installation Instructions: EA01104D →  3

DE Einbauanleitung: EA01104D →  19

Adapter RSE8, RSE5, RSE4

Proline Dosimass, Dosimag

Table of contents

1	Overview of accessories	4
2	Authorized installation personnel	4
3	Safety instructions	4
4	Symbols used	6
5	Dosimass	6
6	Dosimag	12

1 Overview of accessories

The Installation Instructions apply to the following accessories:

Order number	Device component	Contents
50107167	Set, adapter RSE5 8-pin → 4-pin, 4 wire	Adapter for M12×1 socket RSE4/4-pin
50107168	Set, adapter RSE5 8-pin → 5-pin, 5 wire	Adapter for M12×1-socket RSE5/5-pin
50107169	Set, adapter RSE8 8-pin → 8-pin, 5 wire	Adapter for M12×1-socket RSE8/8-pin

2 Authorized installation personnel

Authorization to carry out installation depends on the measuring device's approval type. The table below shows the authorized group of people in each case.

 Whoever carries out the installation has full responsibility to ensure that work is carried out safely and to the required quality standard. He/she must also guarantee the safety of the device following installation.

Measuring device approval	Personnel authorized to carry out installation ¹⁾
No approval	1, 2, 3
With approval (e.g. IECEx)	1, 2, 3
For custody transfer	4

- 1) 1 = Qualified specialist on customer side, 2 = Service technician authorized by Endress+Hauser, 3 = Endress+Hauser (return measuring device to manufacturer)
4 = Check with local approval center if installation/modification must be performed under supervision.

3 Safety instructions

- Only use original parts from Endress+Hauser.
- Comply with national regulations governing mounting, electrical installation, commissioning, maintenance and repair procedures.
- The following requirements must be met with regard to specialized technical staff for the mounting, electrical installation, commissioning, maintenance and repair of the measuring devices:
 - Specialized technical staff must be trained in instrument safety.
 - They must be familiar with the individual operating conditions of the devices.
 - In the case of Ex-certified measuring devices, they must also be trained in explosion protection.
- The measuring device is energized! Risk of fatal injury from electric shock. Open the measuring device only when the device is deenergized.





- When using Ex-certified measuring devices: open device only when in a deenergized state (allow 10 minutes to elapse after switching off the power supply) or in environments which do not have a potentially explosive atmosphere.
- When using measuring devices in safety-related applications in accordance with IEC 61508 or IEC 61511: following installation, re-commission the device in accordance with the Operating Instructions. Document the installation.
- Before removing the device: set the process to a safe state and purge the pipe of dangerous process substances.
- Hot surfaces! Risk of injury! Before commencing work: allow the system and measuring device to cool down to a touchable temperature.
- In the case of devices in custody transfer, the custody transfer status no longer applies once the seal has been removed.
- The Operating Instructions for the device must be followed.
- Risk of damaging the electronic components! Ensure you have a working environment protected from electrostatic discharge.
- After removal of the electronics compartment cover: risk of electrical shock due to missing touch protection!
Turn instrument off before removing internal covers.
- Modifications to the measuring device are not permitted.
- Only open housing for a brief period. Avoid the penetration of foreign bodies, moisture or contaminants.
- Replace defective seals only with original seals from Endress+Hauser.
- If threads are damaged or defective, the measuring device must be repaired.
- Threads (e.g. of the electronics compartment cover and connection compartment cover) must be lubricated if an abrasion-proof dry lubricant is not available. Use acid-free, non-hardening lubricant.
- If, during installation, spacing is reduced or the dielectric strength of the measuring device cannot be guaranteed, perform a test on completion of the work (e.g. high-voltage test in accordance with the manufacturer's instructions).
- Service plug:
 - Do not connect in explosive atmospheres.
 - Only connect to Endress+Hauser service devices.
- Observe the instructions for transporting and returning the device outlined in the Operating Instructions.



If you have any questions, please contact your [Endress+Hauser service organization](#).


4 Symbols used


4.1 Symbols for certain types of information

Symbol	Meaning
	Permitted Procedures, processes or actions that are permitted.
	Forbidden Procedures, processes or actions that are forbidden.
	Tip Indicates additional information.
	Series of steps


5 Dosimass

5.1 Connection with adapter 8 → 8-pin (power supply, pulse output, status output)



 1 Wiring diagram for adapter 8 → 8-pin

- A Socket on device
- B Cable connector
- C Adapter
- 1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 5 (+), pulse output, status output (max. 30 V)
- 6 (-), pulse output (max. 25 mA)
- 7 (-), status output (max. 25 mA)
- 2 Not assigned
- 3 Not assigned
- 8 Not assigned



2 Connection example for adapter 8 → 8-pin (adapter RSE8, 50107169)

A Socket on device

B Cable connector


C Adapter

D PELV or SELV power supply

1 to Pin contacts in plug-in connection


8

5.2 Connection with adapter 8 → 5-pin (power supply, pulse output, status output)



3 Wiring diagram for adapter 8 → 5-pin

- A Socket on device
- B Cable connector
- C Adapter
- 1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 5 (+), pulse output, status output (max. 30 V)
- 6 (-), pulse output (max. 25 mA)
- 7 (-), status output (max. 25 mA)



▣ 4 Connection example for adapter 8 → 5-pin (adapter RSE8, 50107168)

A Socket on device

B Cable connector


C Adapter

D PELV or SELV power supply

1 to Pin contacts in plug-in connection


8

5.3 Connection with adapter 8 → 4-pin (power supply, pulse output, status output)



5 Wiring diagram for adapter 8 → 4-pin

- A Socket on device
- B Cable connector
- C Adapter
- 1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 5 (+), pulse output, status output (max. 30 V)
- 6 (-), pulse output (max. 25 mA)



▣ 6 Connection example for adapter 8 → 4-pin (adapter RSE8, 50107167)

A Socket on device

B Cable connector

C Adapter


D PELV or SELV power supply

1 to Pin contacts in plug-in connection

8

6 Dosimag

6.1 Connection with adapter 8 → 8-pin (power supply, pulse output, status output)



7 Wiring diagram for adapter 8 → 8-pin

A Socket on device

B Cable connector

C Adapter

1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)

4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)

5 (+), pulse output, status output (max. 30 V)


6 (-), pulse output (max. 25 mA)

7 (-), status output (max. 25 mA)

2 Not assigned

3 Not assigned

8 Not assigned



8 Connection example for adapter 8 → 8-pin (adapter RSE8, 50107169)

A Socket on device

B Cable connector


C Adapter

D PELV or SELV power supply

1 to Pin contacts in plug-in connection


8

6.2 Connection with adapter 8 → 5-pin (power supply, pulse output, status output)



9 Wiring diagram for adapter 8 → 5-pin


- A Socket on device
- B Cable connector
- C Adapter
- 1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 5 (+), pulse output, status output (max. 30 V)
- 6 (-), pulse output (max. 25 mA)
- 7 (-), status output (max. 25 mA)



10 Connection example for adapter 8 → 5-pin (adapter RSE8, 50107168)


- A Socket on device
- B Cable connector
- C Adapter
- D PELV or SELV power supply
- 1 to Pin contacts in plug-in connection
- 8

6.3 Connection with adapter 8 → 4-pin (power supply, pulse output, status output)



11 Wiring diagram for adapter 8 → 4-pin

- A Socket on device
- B Cable connector
- C Adapter
- 1 (+), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 4 (-), power supply (24 V DC nominal voltage (20 to 30 VDC), 4.3 W)
- 5 (+), pulse output, status output (max. 30 V)
- 6 (-), pulse output (max. 25 mA)



12 Connection example for adapter 8 → 4-pin (adapter RSE8, 50107167)

- A Socket on device
- B Cable connector
- C Adapter
- D PELV or SELV power supply
- 1 to Pin contacts in plug-in connection
- 8

Adapter RSE8, RSE5, RSE4

Proline Dosimass, Dosimag

Inhaltsverzeichnis

1	Übersicht Zubehörteile	20
2	Einbauberechtigte Personen	20
3	Sicherheitshinweise	20
4	Verwendete Symbole	22
5	Dosimass	22
6	Dosimag	28


1 Übersicht Zubehörteile

Die Einbauanleitung ist für folgende Zubehörteile gültig:

Bestellnummer	Gerätekomponente	Inhalt
50107167	Set Adapter RSE5 8pol → 4pol, 4 Leiter	Adapter auf M12×1-Buchse RSE4/4pol
50107168	Set Adapter RSE5 8pol → 5pol, 5 Leiter	Adapter auf M12×1-Buchse RSE5/5pol
50107169	Set Adapter RSE8 8pol → 8pol, 5 Leiter	Adapter auf M12×1-Buchse RSE8/8pol

2 Einbauberechtigte Personen

Die Berechtigung zur Durchführung eines Einbaus ist von der Zulassung des Messgeräts abhängig. Die Tabelle zeigt den jeweils berechtigten Personenkreis.

 Die Person, die einen Einbau vornimmt, übernimmt die Verantwortung für die Sicherheit während der Arbeiten, die Qualität der Ausführung und die Sicherheit des Geräts nach dem Einbau.

Zulassung des Messgeräts	Einbauberechtigter Personenkreis ¹⁾
Ohne Zulassung	1, 2, 3
Mit Zulassung (z.B. IECEx)	1, 2, 3
Bei eichfähigem Verkehr	4

- 1) 1 = Ausgebildete Fachkraft des Kunden, 2 = Von Endress+Hauser autorisierter Servicetechniker,
 3 = Endress+Hauser (Messgerät an Hersteller zurücksenden)
 4 = Mit der lokalen Zulassungsstelle prüfen, ob ein Ein-/Umbau unter Aufsicht erfolgen muss.

3 Sicherheitshinweise

- Nur Originalteile von Endress+Hauser verwenden.
- Nationale Vorschriften bezüglich der Montage, elektrischen Installation, Inbetriebnahme, Wartung und Reparatur einhalten.
- Folgende Anforderungen an das Fachpersonal für Montage, elektrische Installation, Inbetriebnahme, Wartung und Reparatur der Messgeräte müssen erfüllt sein:
 - In Gerätesicherheit ausgebildet.
 - Mit den jeweiligen Einsatzbedingungen der Geräte vertraut.
 - Bei Ex-zertifizierten Messgeräten: zusätzlich im Explosionsschutz ausgebildet.
- Messgerät unter Spannung! Lebensgefahr durch Stromschlag. Messgerät nur im spannungslosen Zustand öffnen.
- Bei Ex-zertifizierten Messgeräten: Nur in spannungslosem Zustand (nach Berücksichtigung einer Wartezeit von 10 Minuten nach Abschalten der Energiezufuhr) oder in Umgebungen öffnen, die keine explosionsfähige Atmosphäre enthalten.





- Bei Messgeräten in sicherheitstechnischen Applikationen gemäß IEC 61508 bzw. IEC 61511: Nach Einbau Neuinbetriebnahme gemäß Betriebsanleitung durchführen. Einbau dokumentieren.
- Vor einem Geräteausbau: Prozess in sicheren Zustand bringen und Leitung von gefährlichen Prozessstoffen befreien.
- Verbrennungsgefahr durch heiße Oberflächen! Vor Arbeitsbeginn: Anlage und Messgerät auf berührungssichere Temperatur abkühlen.
- Bei Messgeräten im abrechnungspflichtigen Verkehr: Nach Entfernen der Plombe ist der geeichte Zustand aufgehoben.
- Die Betriebsanleitung zum Messgerät ist zu beachten.
- Beschädigungsgefahr elektronischer Bauteile! Eine ESD-geschützte Arbeitsumgebung herstellen.
- Nach Entfernen der Elektronikabdeckung: Stromschlaggefahr durch aufgehobenen Berührungsschutz!
Messgerät ausschalten, bevor interne Abdeckungen entfernt werden.
- Änderungen am Messgerät sind nicht zulässig.
- Gehäuse nur kurzzeitig öffnen. Eindringen von Fremdkörpern, Feuchtigkeit oder Verunreinigung vermeiden.
- Defekte Dichtungen nur durch Original-Dichtungen von Endress+Hauser ersetzen.
- Defekte Gewinde erfordern eine Instandsetzung des Messgeräts.
- Gewinde (z.B. von Elektronikraum- und Anschlussraumdeckel) müssen geschmiert sein, sofern keine abriebfeste Trockenschmierung vorhanden ist. Säurefreies, nicht härtendes Fett verwenden.
- Wenn bei den Einbauarbeiten Abstände reduziert oder die Spannungsfestigkeit des Messgeräts nicht sichergestellt werden kann: Prüfung nach Abschluss der Arbeiten durchführen (z.B. Hochspannungstest gemäß Herstellerangaben).
- Servicestecker:
 - Nicht in explosionsfähiger Atmosphäre anschließen.
 - Nur an Servicegeräte von Endress+Hauser anschließen.
- Die in der Betriebsanleitung aufgeführten Hinweise zum Transport und zur Rücksendung beachten.



Bei Fragen kontaktieren Sie bitte Ihre zuständige [Endress+Hauser Serviceorganisation](#).


4 Verwendete Symbole


4.1 Symbole für Informationstypen

Symbol	Bedeutung
	Erlaubt Abläufe, Prozesse oder Handlungen, die erlaubt sind.
	Verboten Abläufe, Prozesse oder Handlungen, die verboten sind.
	Tip Kennzeichnet zusätzliche Informationen.
	Handlungsschritte


5 Dosimass

5.1 Anschluss mit Adapter 8 → 8pol (Energieversorgung, Impuls-, Statusausgang)



 13 Anschlussplan Adapter 8 → 8pol


- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- 1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 5 (+), Impuls-, Statusausgang (max. 30 V)
- 6 (-), Impulsausgang (max. 25 mA)
- 7 (-), Statusausgang (max. 25 mA)
- 2 nicht belegt
- 3 nicht belegt
- 8 nicht belegt



14 Anschlussbeispiel Adapter 8 → 8pol (Adapter RSE8, 50107169)

- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung

5.2 Anschluss mit Adapter 8 → 5pol (Energieversorgung, Impuls-, Statusausgang)



15 Anschlussplan Adapter 8 → 5pol

A Steckbuchse am Gerät

B Kabelstecker

C Adapter


1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)

4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)

5 (+), Impuls-, Statusausgang (max. 30 V)

6 (-), Impulsausgang (max. 25 mA)


7 (-), Statusausgang (max. 25 mA)



▣ 16 Anschlussbeispiel Adapter 8 → 5pol (Adapter RSE8, 50107168)


- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung

5.3 Anschluss mit Adapter 8 → 4pol (Energieversorgung, Impuls-, Statusausgang)



17 Anschlussplan Adapter 8 → 4pol

- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- 1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 5 (+), Impuls-, Statusausgang (max. 30 V)
- 6 (-), Impulsausgang (max. 25 mA)




▣ 18 Anschlussbeispiel Adapter 8 → 4pol (Adapter RSE8, 50107167)

- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung

6 Dosimag

6.1 Anschluss mit Adapter 8 → 8pol (Energieversorgung, Impuls-, Statusausgang)



▣ 19 Anschlussplan Adapter 8 → 8pol

A Steckbuchse am Gerät

B Kabelstecker

C Adapter

1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)

4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)

5 (+), Impuls-, Statusausgang (max. 30 V)


6 (-), Impulsausgang (max. 25 mA)

7 (-), Statusausgang (max. 25 mA)

2 nicht belegt

3 nicht belegt


8 nicht belegt



20 Anschlussbeispiel Adapter 8 → 8pol (Adapter RSEB, 50107169)


- A Steckbuchse am Gerät
 - B Kabelstecker
 - C Adapter
 - D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung

6.2 Anschluss mit Adapter 8 → 5pol (Energieversorgung, Impuls-, Statusausgang)



21 Anschlussplan Adapter 8 → 5pol


- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- 1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 5 (+), Impuls-, Statusausgang (max. 30 V)
- 6 (-), Impulsausgang (max. 25 mA)
- 7 (-), Statusausgang (max. 25 mA)



▣ 22 Anschlussbeispiel Adapter 8 → 5pol (Adapter RSE8, 50107168)


- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung

6.3 Anschluss mit Adapter 8 → 4pol (Energieversorgung, Impuls-, Statusausgang)



▣ 23 Anschlussplan Adapter 8 → 4pol

- A Steckbuchse am Gerät
- B Kabelstecker
- C Adapter
- 1 (+), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 4 (-), Energieversorgung (24 V DC Nominalspannung (20...30 VDC), 4,3 W)
- 5 (+), Impuls-, Statusausgang (max. 30 V)
- 6 (-), Impulsausgang (max. 25 mA)



▣ 24 Anschlussbeispiel Adapter 8 → 4pol (Adapter RSE8, 50107167)

- A Steckbuchse am Gerät
 - B Kabelstecker
 - C Adapter
 - D PELV oder SELV Spannungsversorgung
- 1...8 Kontaktstifte in der Steckverbindung



71409289

www.addresses.endress.com
