

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser

Obere Wank 1, 87484 Nesselwang, Germany

the Certificate No: **Z01604** for the PROFIBUS Device:

Model Name: iTEMP TMT162 3.02
Revision: 02; SW/FW: 01.01.03; HW: 01.00.07
GSD: EH021549.GSD; File Version 2.0
PA139700.gsd; PA139701.gsd
PA139702.gsd; PA139703.gsd

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	DP-V0	MS0, Freeze, Set_Slave_Add
<input checked="" type="checkbox"/>	DP-V1	MS1, MS2, I&M
<input checked="" type="checkbox"/>	Profile	PA Devices Profile 3.02
<input checked="" type="checkbox"/>	Physical Layer	MBP

Test Report Number: itm 3007 PA 01
Authorized Test Laboratory: itm, München, Germany

The tests were executed in accordance with the following documents:
"Test Specifications for PROFIBUS DP Slaves, Version 3.0 from November 2005" and
"Test Specifications for PA Devices, Profile 3.02, Version 5.0.4, February 2010".

This certificate is granted according to the document:
"Framework for testing and certification of PROFIBUS and PROFINET products".


For all products that are placed in circulation by **April 11, 2023** the certificate is valid for life.

Karlsruhe, April 03, 2020

Board of PROFIBUS Nutzerorganisation e. V.




(Official in Charge)



(Karsten Schneider)





(Dr. Jörg Hähnliche)