

# **Certificate of Compliance**

Certificate:

1843172 (LR 53988)

**Master Contract:** 

151079

Project:

1843172

Date Issued:

2007/01/19

Issued to:

Endress + Hauser GmbH + Co. KG

Haupstrasse 1

Maulburg, Baden-Wurttemberg 79690

Germany

Attention: Stefan Muth

The products listed below are eligible to bear the CSA Mark shown



Issued by:

Edward Foo, C.E.T.

Authorized by: Patricia Pasemko, Operations

Manager

tatinia Dasent)

#### **PRODUCTS**

CLASS 2252 05 - PROCESS CONTROL EQUIPMENT

• Pressure transducer, stationary mounted, Model Cerabar T PMC 131-C\*\*\*, rated 11-30Vdc, 4-20mA 2-wire output, Ta = -20 to +85°C. Encl. Type 4X (plug connector) or 6P (cable version). Pollution Degree 2; Installation category II environment.

#### Notes:

- The model suffixes asterisks "\*\*\*" denote alpha-numeric characters representing model variations;
- Device powered by SELV or Class 2 supply source;



Certificate:

1843172 (LR 53988)

**Master Contract:** 

151079

Project:

1843172

**Date Issued:** 

2007/01/19

- Device for use in equipment where the suitability of the combination is to be determined by Canadian Standards Association International.

### **APPLICABLE REQUIREMENTS**

CAN/CSA-C22.2 No. 61010-1-04 - Safety Requirements for Electrical Equipment, for Measurement, Control, and Laboratory Use, Second Edition, Part 1: General Requirements.

DQD 507 Rev. 2004-06-30



# Supplement to Certificate of Compliance

Certificate:

1843172

**Master Contract:** 

151079

The products listed, including the latest revision described below, are eligible to be marked in accordance with the referenced Certificate.

## **Product Certification History**

Project	Date	Description
1843172	2007/01/19	CSA Certification of Pressure Transducer for Ordinary Locations