


Device Registration

Manufacturer	Endress+Hauser Flowtec AG Kaegenstrasse 7 Reinach, BL, 4153 Switzerland	
Model Type	Prowirl 200 2-wire Vortex Flowmeter	
H1 ITK Profile	6	
Manufacturer Identification (MANUFAC_ID)	452B48	
Device Type (DEV_TYPE)	1038	
Device Revision (DEV_REV)	01	
Stack Communication Profiles	31PS, 32LT	
H1 Physical Layer Profiles	111, 511	
Device Test Campaign	IT094200	
Stack Test Campaign	CT0177FF	
H1 Physical Layer Test Report	PT-445	
Device Descriptions	0101.FFO 0101.SYM 0101.FF5 0101.SY5 010101.CFF	(CRC: 0x6350B7FB) (CRC: 0x4B46410C) (CRC: 0x4775573E) (CRC: 0x4543D8D5) (CRC: 0xC9BFD6B5)
Capability File		
Mandatory Features	Resource Block Alarms and Events Function Block Linking Trending Multi-Bit Alert Reporting Field Diagnostics	
Function Blocks	Analog Input Discrete Input Proportional Integral Derivative Integrator Multiple Analog Output Multiple Discrete Output	
Registration Number	IT/094200/1	
Registration Program	FF-524 2.0	
Date Issued	2014-03-26	

**Authorized
Mark**

Use of mark subject to
Graphic Identify Standards Guide
(MT-042)

Richard J. Timoney
President and Chief Executive Officer