



Device Registration

Manufacturer	Endress+Hauser Flowtec AG Kaegenstrasse 7 Reinach, BL, 4153 Switzerland	
Model Type	Promag 200 2-wire Magnetic Flowmeter	
H1 ITK Profile	6	
Manufacturer Identification (MANUFAC_ID)	0x452B48	
Device Type (DEV_TYPE)	0x1048	
Device Revision (DEV_REV)	0x01	
Stack Communication Profiles	31PS, 32LT	
H1 Physical Layer Profiles	111, 511	
Device Test Campaign	IT102900	
Stack Test Campaign	CT0175FF	
H1 Physical Layer Test Report	PT-492	
Device Descriptions	0101.FFO	(CRC: 0x179710D)
	0101.SYM	(CRC: 0x8101EFF4)
	0101.FF5	(CRC: 0x181DCC20)
	0101.SY5	(CRC: 0x191D8382)
Capability File	010101.CFF	(CRC: 0x9263B29B)
Mandatory Features	Resource Block Alarms and Events Function Block Linking Trending Multi-Bit Alert Reporting Field Diagnostics	
Function Blocks	Analog Input Discrete Input Proportional Integral Derivative Integrator Multiple Discrete Output	
Registration Number	IT/102900/1	
Registration Program	FF-524 2.11	
Date Issued	2015-10-27	

**Authorized
Mark**

Use of mark subject to
Graphic Identify Style Guide



T. J. Masters

Ted Masters

President and Chief Executive Officer