

TM



FOUNDATION

FIELDBUS FOUNDATION

DD/CF

REGISTRATION

Presented To: Endress+Hauser Process Solutions AG

Name: Flowmeter

Model: Prowin73

Device Type: 0x1057

Device Revision: 0x01

Original Interoperability Test Campaign: IT022900

Device Description: 0101.FFO (CRC32:4314C8ED)

: 0101.SYM (CRC32:FB95883)


Capabilities File: 010102.CFF (CRC32:3BB26DF2)

ITK Version: 4.60

DD/CF Test Campaign Number: IT029700

DD/CF Re-Registration Date: 11/16/2004

The above device has successfully completed rigorous testing by the Fieldbus Foundation and has received registration and the right to use the FF checkmark logo as specified by MT-045.



Heather Cannard

Heather Cannard
Test Technician

Richard J. Timoney

Richard J. Timoney
President