

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser Wetzer GmbH & Co. KG

Obere Wank 1, 87484 Nesselwang, Germany

the Certificate No: **Z60003** for the PROFINET PA Profile Device:

Model Name:

iTEMP TMT86

Revision:

SW/FW: V1.00.36; HW: 100

Device ID:

0xB300

GSD:

GSDML-V2.43-PA_Profile_V4.02-B300-TEMPERATURE-20230217.XML

DAP:

ID DAP_APL_10MBITBP, 0x00FD0000

This certificate confirms that the product has successfully passed the certification tests with the following scope:

☑ PA Profile	V4.02MU1
☐ PA_Tester_Version	V2.43.3.3
☑ Tester	ifak e.V., Magdeburg, Germany; PA314-01

This certificate is granted according to the document:

"Framework for testing and certification of PROFIBUS and PROFINET products".

For all products that are placed in circulation by **November 07, 2026** the certificate is valid for life.

Karlsruhe, November 09, 2023

(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.

ave filish

(Xaver Schmidt)

(Frank Moritz)

