

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser Flowtec AG

Christoph Merian-Ring 4, 4153 Reinach, BL, Switzerland

the Certificate No: **Z60016** for the PROFINET PA Profile Device:

Model Name:

Prowirl 200 PROFINET-APL

Revision:

SW/FW: V1.0.2; HW: 103

Device ID:

0xB330

GSD:

GSDML-V2.43-EH-Prowirl_200_APL-20240409.xml

DAP:

DAP 1 V1.00.01, 0x00fd0000

This certificate confirms that the product has successfully passed the certification tests with the following scope:

☑ PA Profile	V4.02MU1
☐ PA_Tester_Version	Version 2.44.4
☑ Tester	ifak, Magdeburg, Germany; PA326-01

This certificate is granted according to the document:

"Framework for testing and certification of PROFIBUS and PROFINET products".

For all products that are placed in circulation by June 17, 2027 the certificate is valid for life.

Karlsruhe, October 21, 2024

(Official in Charge)

PROFIprocess

Board of PROFIBUS Nutzerorganisation e. V.

ave filish

(Xaver Schmidt)

(Frank Moritz)