

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser Flowtec AG

Christoph Merian-Ring 4, 4153 Reinach, BL, Switzerland

the Certificate No: **Z60015** for the PROFINET PA Profile Device:

Model Name: Prowirl 200 PROFINET-APL
Revision: SW/FW: V1.0.2; HW: 103
Device ID: 0xB330
GSD: GSDML-V2.43-PA_Profile_V4.02-B330-FLOW_VORTEX-20230217.xml
DAP: ID_DAP_APL_10MBITBP, 0x00FD0002

This certificate confirms that the product has successfully passed the certification tests with the following scope:

☒ PA Profile	V4.02MU1
☒ PA_Tester_Version	Version 2.44.4
☒ Tester	ifak, Magdeburg, Germany; PA326-01


This certificate is granted according to the document:

"Framework for testing and certification of PROFIBUS and PROFINET products".

For all products that are placed in circulation by **June 17, 2027** the certificate is valid for life.

Karlsruhe, October 21, 2024

Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)



(Xaver Schmidt)



(Frank Moritz)

