



EL Class I

Date of issue: 3 May 2021

Valid until: 31 December 2025

EHEDG hereby declares that the product


Liquipoint (O)FTW23 with PEEK sensor tip and EPDM gasket


from

Endress+Hauser SE+Co.KG, Hauptstrasse 1 , 79689 Maulburg, Germany

*has/have been evaluated for compliance and meets/meet the current criteria for
Hygienic Equipment Design of the EHEDG*

Certificate No. EHEDG-C2100020

Signed  President EHEDG
Hein Timmerman

Signed  EHEDG Certification Officer
Karlijn Faber

EHEDG
Karspeldreef 8
1101 CJ Amsterdam
Netherlands

©EHEDG