

## Manufacturer Declaration

Endress+Hauser GmbH+Co. KG, Hauptstraße 1, D-79689 Maulburg

Being the manufacturer, declares that the product stated below

**Liquipoint FTW23 with  
Option "7"- IO-Link; DC-PNP**

is in conformity with the following provisions:



IO-Link Interface and System Specification V1.1

Device Test System Revision: 3.0.9.31; Test ID: 562067926

The conformity tests are documented in test report:


IOL-Test\_FW\_02.00.00-038\_2017-10-23.pdf

*This declaration of compliance is exclusively valid for the listed products and accessory in delivery status. The validity of this document for placing the product on the market expires 3 years after the date of issue.*

Maulburg, 7-November-2017  
Endress+Hauser GmbH+Co. KG



i. A. Armin Hummelbrumm  
Group Leader R&D  
Dept. TPS



i. V. Dietmar Frühauf  
Dept. Manager Level Switch  
Research & Development Devices