

Manufacturer Declaration

Endress+Hauser GmbH+Co. KG, Hauptstraße 1, D-79689 Maulburg

Being the manufacturer, declares that the product stated below

**Cerabar PMP23 with
Option "7"- IO-Link; 4-20mA**

is in conformity with the following provisions:



IO-Link Interface and System Specification V1.1

Device Test System Revision: 3.0.9.31; Test ID: 562080903


The conformity tests are documented in test report:

IOL-Test_FW_01.00.02_2017-10-23.pdf

This declaration of compliance is exclusively valid for the listed products and accessory in delivery status. The validity of this document for placing the product on the market expires 3 years after the date of issue.

Maulburg, 7-November-2017
Endress+Hauser GmbH+Co. KG


i. A. Armin Hummelbrumm
Group Leader R&D
Dept. TPS


i. V. Marc Schlachter
Dept. Manager Pressure
Research & Development Devices