

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser

Obere Wank 1, 87484 Nesselwang, Germany

the Certificate No: **Z11306** for the PROFINET Device:

Model Name: Memograph M RSG45
 Revision: SW/FW: V2.04.02; HW: V1
 Identnumber: 0x0011; 0x86FA
 GSD: GSDML-V2.34-EH-RSG45-20180522.xml
 DAP: DAP, 0x80010000

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	PNIO_Version	V2.34
<input checked="" type="checkbox"/>	Conformance Class	B
<input checked="" type="checkbox"/>	Optional Features	Legacy, MRP, DeviceAccess
<input checked="" type="checkbox"/>	Netload Class	III
<input checked="" type="checkbox"/>	PNIO_Tester_Version	Version 2.34.0
<input checked="" type="checkbox"/>	Tester	AIT, Gummersbach, Germany; PN AIT-18-07


This certificate is granted according to the document:

“Framework for testing and certification of PROFIBUS and PROFINET products”.

For all products that are placed in circulation by **June 26, 2021** the certificate is valid for life.

Karlsruhe, September 20, 2018

Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)




(Karsten Schneider)



(Dr. Jörg Hähnliche)