

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser Conducta GmbH+Co. KG
Dieselstrasse 24, 70839 Gerlingen, Germany

the Certificate No: **Z12698** for the PROFINET Device:

Model Name: Liquiline CM44x
Revision: SW/FW: V1.8.0; HW: 202
Identnumber: 0x0011; 0x859C
GSD: GSDML-V2.34-EH-CM44x-20200717.xml
DAP: DAP 1, 0x00FD0000

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/> PNIO_Version	V2.34
<input checked="" type="checkbox"/> Conformance Class	B
<input checked="" type="checkbox"/> Optional Features	Legacy, DeviceAccessSupported
<input checked="" type="checkbox"/> Netload Class	II
<input checked="" type="checkbox"/> PNIO_Tester_Version	Version 2.4.1
<input checked="" type="checkbox"/> Tester	AIT, Gummersbach, Germany; PN AIT-20-14

This certificate is granted according to the document:

“Framework for testing and certification of PROFIBUS and PROFINET products”.

For all products that are placed in circulation by **August 05, 2026** the certificate is valid for life.

Karlsruhe, April 22, 2024

Board of PROFIBUS Nutzerorganisation e. V.



(Official in Charge)



(Xaver Schmidt)




(Frank Moritz)