

Company **Endress+Hauser Flowtec AG, Kägenstrasse 7, CH-4153 Reinach**
declares as manufacturer, for the following product


Product **Proline Promag 100, Proline Promass 100, Proline Prosonic Flow 100**
Proline Promag 200, Proline Promass 200, Proline Prosonic Flow 200, Proline Prowirl 200
Proline Promag 300, Proline Promass 300
Proline Promag 400
Proline Promag 500, Proline Promass 500
with Heartbeat Technology

Statement All the above-mentioned devices with Heartbeat Technology manufactured by Endress+Hauser Flowtec AG which, according to the attestation by TÜV SÜD Industrie Service GmbH, meet the requirements for traceable verification according to DIN EN ISO 9001:2008 - Section "7.6 a) Control of monitoring and measuring equipment" also meet the requirement for traceable verification according to ISO 9001:2015, Section "7.1.5 Monitoring and measuring resources" and Section "7.1.5.2 Measurement Traceability".

Attestations for Proline devices with Heartbeat Technology issued by TÜV SÜD Industrie Service GmbH and referencing DIN EN ISO 9001:2008 as the Test specification remain valid without restrictions.

This declaration applies only for the listed product(s) in delivery condition.

Reinach, 30. April 2024
Endress+Hauser Flowtec AG



Dr. Mirko Lehmann
Managing Director
Endress+Hauser Flowtec AG



Dipl.-Ing. Michael Karolzak
Principal Expert Functional Safety
Endress+Hauser Flowtec AG