



Certificate of Registration HCF Verified

<u>Endress + Hauser</u> Manufacturer	<u>t-mass L/T 150</u> Product Name
<u>0011</u> Manufacturer ID (Hex)	<u>1168</u> Expanded Device Type (Hex)
<u>6</u> HART Protocol Revision	<u>01</u> Device Revision (Hex)
<u>01</u> Hardware Revision (Hex)	<u>01</u> Software Revision (Hex)
<u>6/12/2014</u> Test Date	<u>HCF</u> Verification Method

The above product has successfully completed the validation process and meets the requirements to be "HART REGISTERED".

"HART REGISTERED" products conform to GB/T 29910.1-6-2013 and IEC 61158 standards.

Registration
Number: L2-06-1000-337

Registration
Issue Date: June 12, 2014

HCF QA
Approval: *J. F. Masters*

