

CERTIFICATE OF COMPLIANCE



*TNO Certification
hereby declares that the product*

Liquiphant T, type T

From

Endress+Hauser GmbH+Co. KG

*has been evaluated for compliance with the Hygienic Equipment Design Criteria
of the EHEDG, Document No. 8, by:*

*TNO Nutrition and Food Research at Zeist, Netherlands
and meets the criteria of this document as demonstrated by:*

Evaluation Report No. V5496

Signed

Evaluation Officer

Date December 07, 2003

Signed

Managing Director, TNO Certification

Date December 07, 2003



Certificate No. C03-5617

*TNO Certification BV, P.O. Box 541, 7300 AM Apeldoorn, Netherlands
©EHEDG*



-Registration Liquiphant T FTL20H

Certificate No.: C03-5617 / Report No. V5496

Overview registered process connections

Type	Process connection	Order code: 020
FTL20H	E+H flush mounted for 1" welding neck	UPJ including accessory weld-in adapter 1" Part No. 52001047 or 52006909
	Thread ISO228 G 3/4 G1	GDJ GEJ including accessory weld-in adapter 3/4" Part No. 52018765 1" Part No. 52001051 or 52011896
	Clamp ISO2852 2"	TDJ