



Certificate

PROFIBUS Nutzerorganisation e.V. grants to

Endress+Hauser GmbH & Co. KG
Hauptsraße 1; D-79689 Maulburg

the Certificate No.: **Z01251** for the PROFIBUS Slave:

Product Name: Prosonic S FMU95
Revision: 01.00; SW/FW: All; HW: All
GSD: EH3x154E.gsd; File Version: 0.0 2006/08/09
PA03970F.gsd

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/>	DP-V0	MS0, Freeze, Set_Slave_Add
<input checked="" type="checkbox"/>	DP-V1	MS2, I&M
<input type="checkbox"/>	DP-V2	
<input checked="" type="checkbox"/>	Profile	PA Devices V 3.0
<input checked="" type="checkbox"/>	Physical Layer	RS485

Test Report Number: itm 648 PA 01/01
Authorized Test Laboratory: itm, München, Germany
Expiry date of Certificate: March 08, 2010


The tests were executed in accordance with the following documents:
"Test Specifications for PROFIBUS DP Slaves, Version 3.0 from November 2005" and
"Test Specification for PA Devices, Profile 3.0; Version 4.1, November 2004".
This certificate is granted according to the document "Framework for testing and certification of PROFIBUS products".

Karlsruhe, May 04, 2007





.....
(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.



(E. Küster)



(Prof. K. Bender)