

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

**Endress+Hauser Flowtec AG,
Kägenstraße 7, 4153 Reinach, Switzerland**

the Certificate No: **Z00869** for the PROFIBUS Slave:

Model Name: PROSONIC FLOW 93 PA
Revision: Profile 3.0; SW/FW: All 2; HW: All 2
GSD: eh3x1530.GSD; File Version: 0.3
 PA139741.GSD

This certificate confirms that the product has successfully passed the certification tests with the following scope:


<input checked="" type="checkbox"/>	DP-V0	MS0, Sync, Freeze, Set_Slave_Add
<input checked="" type="checkbox"/>	DP-V1	MS1, MS2, I&M
<input checked="" type="checkbox"/>	Profile	PA Devices V3.0
<input checked="" type="checkbox"/>	Physical Layer	MBP

Test Report Number: 673 PA 01/03
Authorized Test Laboratory: itm, München, Germany

The tests were executed in accordance with the following documents:
 "Test Specifications for PROFIBUS DP Slaves, Version 3.0 from November 2005",
 "Test Specification for PA Devices, Profile 3.0, Version 4.1 from November 2004.

This certificate is granted according to the document:
 "Framework for testing and certification of PROFIBUS products".

For all products that are placed in circulation by November 11, 2012 the certificate is valid for life.




(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.



(Jörg Freitag)

(K.-P. Lindner)