


KA 079R/09/c4/07.98
No.: 50092424




Omnigrad TMT147-TMT146

2 Draht Messumformer für
Thermowiderstände Pt100
2 Wire transmitter for
thermoresistances Pt100
Trasmettitore a 2 fili per
termoresistenze Pt100
Trasmetteur liaison 2 fils pour
thermoresistances Pt100

Betriebsanleitung
Installation Manual
Manuale di installazione
Mise en service



Endress+Hauser 



TMT147-R TMT146-R

- A** Pt100 2-Draht
Pt100 2-wires
Pt100 2-fili
Pt100 2 fils
- B** Pt100 3-Draht
Pt100 3-wires
Pt100 3-fili
Pt100 3 fils
- C** Messspanne
Measuring range
Campo di misura
Etendue de mesure

DS3P013A



Endress+Hauser

